

Characterization of Antibody Drug Conjugates using Microfluidic CE-MS

Erin A. Redman, Joshua P. Guerrette, Michael P. Goodwin, J. Scott Mellors
908 Devices Inc., Boston, MA 02210

Overview


- ▶ A microfluidic CE-MS interface was used for the characterization of biotherapeutic monoclonal antibodies and ADC mimics. Multiple levels of analysis were demonstrated to characterize the molecules.
- ▶ Charge variants of the intact molecules were separated and analyzed via on-line ESI-MS. After conjugation, the complexity of the IgG increased significantly but DAR species and other variants of the ADC mimics were characterized without additional sample prep.
- ▶ Additional middle-up analysis of the IgG and ADC mimics revealed further information about the extent of conjugation and the localization of the modified residues.
- ▶ All types of analyses were complete in less than 5 minutes indicating that the microfluidic interface is a simple, rapid means of characterizing complex mAb based therapeutics.

Introduction

Capillary electrophoresis (CE) has been demonstrated to be effective for analyzing complex biological samples, such as monoclonal antibodies (mAbs) and other protein based therapies. When interfaced with mass spectrometry (MS) analysis the integrated techniques present a powerful platform for characterizing biomolecules. Utilizing microfluidic technology to perform CE-MS analysis enables much faster and more efficient separations than can be achieved using capillary tubes. The work described here utilizes a microchip CE-MS interface for the characterization of IgGs and antibody drug conjugate (ADC) mimics.

Methods

ZipChips™ utilize microfluidic technology to harness the inherent speed and efficiency of zone electrophoresis separations. The device design incorporates an injection cross, serpentine separation channel, and an integrated ESI emitter where electrospray is generated directly off the corner of the device. Highly uniform and stable surface coatings suppress the electroosmotic flow and yield highly efficient separations.



Sample prep. ADC mimics were generated using commercially available antibody dye conjugation kits (Thermo Scientific). For analysis, intact IgGs and ADC mimics were diluted to 0.5 mg/mL with LC/MS water. No further sample prep was necessary.

After dilution to 0.5 mg/mL, reduction of the intact molecules was performed with 10 mM DTT in 50 mM ammonium bicarbonate buffer (pH~7). Samples were incubated at 37C for 30 minutes.




IgGs and ADC mimics were digested with the IdeS protease (Fabricator, Genovis) according to manufacturer specifications in 50 mM ammonium bicarbonate buffer (pH~6.5). Digested molecules were then reduced with 10 mM DTT and diluted to 0.5 mg/mL for analysis.

Analysis. ZipChipHR chips were used for all analyses. Intact analysis was performed using the Intact Antibody Analysis Kit (908 Devices, Inc.). Middle-up analysis was performed using a mixture of methanol/water/formic acid. All separations were performed at 500 V/cm.

Data Collection and Processing. An Exactive Plus EMR was used for the intact mAb and ADC characterization. Middle-up analysis of the mAbs and ADC mimics was performed on a Thermo LTQ-XL. Biopharmafinder 2.0 was used for data processing. Images were generated using Igor Pro 7 (WaveMetrics, Inc.).




Unconjugated IgG

This analysis technique is capable of separating charge variants of intact monoclonal antibodies while also yielding high quality MS data in real time. Interfacing the chip separation with online MS analysis enhances the capabilities of the MS and enables rapid and thorough characterization of biotherapeutic products.






IgG Conjugated with Cationic Drug Mimic

ADC mimics were made using commercially available NHS-Dye reagents to create linkages at solvent accessible lysine residues. The non-specific nature and thoroughness of the conjugation reaction resulted in very complex mimics.



IgG Conjugated with Neutral Drug Mimic

Conjugation with the neutral drug mimic had a greater impact on the electrophoretic mobility of the IgG molecule causing the DAR species to elute over a greater timeframe.



Conclusions

The microfluidic interface is an effective tool for characterizing biotherapeutic proteins. The interface can be used to generate high efficiency separations of intact mAb variants with online MS analysis for thorough characterization of mAb and complex ADC samples. Reduced analysis and digestion with IdeS protease was demonstrated for obtaining chain specific modifications and additional detail about conjugation. The short analysis times for both intact analysis and middle-up analysis make this technique amenable to implementation as a routine method for characterizing molecules and monitoring key product attributes.