

Rapid, Multi-level Analysis of Complex Proteins by Microchip Capillary Electrophoresis-ESI-MS

Adi Kulkarni, Erin A. Redman, Ashley Bell, J. Scott Mellors
908 Devices, Inc., Boston, MA 02210

Overview

Mass spec characterization of proteins often requires multiple “levels” of analysis to generate a complete picture of the molecule. These levels are commonly referred to as top-down, middle-down/up, and bottom-up. While LC-MS methods have been developed for each of these levels, those methods typically require completely different columns and mobile phases for the different levels; and the methods typically yield relatively poor results for middle and top down analysis. Microchip CE-ESI-MS can achieve excellent performance for all of these levels, using the same exact experimental conditions. This enables rapid, back-to-back multilevel characterization of a protein, with no down-time between samples. This study demonstrates that capability for a variety of proteins, including complex glycoproteins and monoclonal antibodies.

Methods

Sample Preparation and Analysis. Proteins were prepared for intact analysis by simply diluting to an appropriate concentration before loading into the microchip. Middle-down analysis was performed on monoclonal antibodies by subjecting them to limited proteolysis using the IDES enzyme and reducing them with DTT. Bottom-up analysis utilized a standard trypsin digestion method in ammonium bicarbonate buffer.



Instrumentation. All work was performed using a commercially available microfluidic CE-ESI system (ZipChip, 908 Devices Inc.). The microfluidic devices utilize a covalently attached, neutral polymer surface coating to prevent analyte interactions and suppress electroosmotic flow. For the work shown here, a chip with a 22 cm long separation channel was used (ZipChip HR). Data were collected on a Thermo Tribrid Fusion Lumos mass spectrometer or a Thermo Exactive Plus EMR mass spectrometer. To achieve successful analysis at all three structural levels, a BGE with a relatively low pH was chosen. The ZipChip “Metabolite BGE” (908 Devices Inc.) is a low pH mixture containing water, methanol, and formic acid. While it was developed for optimizing metabolite separations, it is well suited for many peptide and intact protein analyses as well.

Data Processing. Data were visualized using Thermo Xcalibur QualBrowser. Intact and bottom-up data were processed using Thermo BioPharmaFinder 1.0. Top-down fragmentation spectra were deconvoluted using the Xtract algorithm in Thermo FreeStyle. Residue cleavage was calculated using ProSight Lite 1.4 (Northwestern University).


Device Schematic




Intact Analysis




Bottom-up Analysis



Middle-up Analysis



Top-down Analysis



Conclusions

Much of our previous work has focused on fully optimizing methods to achieve the best possible results for different applications of the microfluidic CE-MS (ZipChip) platform. If we instead consider how much can be accomplished without changing either the chip or the BGE, we see a different advantage of this platform. The ability to run such a wide range of applications without making any changes to the system means that a wide variety of samples can simply and efficiently be run on the mass spec. For labs that commonly deal with a variety of samples, this type of flexibility can save a lot of time and allow maximum use of the mass spec. It's also worth noting that while not shown on this poster, the background electrolyte used for this work was originally optimized for metabolite analysis.

The technologies discussed in this poster are the subject of one or more granted/pending patents. www.908devices.com/patents/