

TECH NOTE 1.0

Next-Generation Handheld High-Pressure Mass Spectrometry (HPMS) with an APCI Dual-Polarity Source for Threat Detection


CLASSES	SUBCLASSES	EXEMPLAR TARGETS	DEVICE CAPABILITY
Chemical Warfare Agents (CWAs)	Nerve	G and V series	BB Demo
	Blister	HD, L	
Toxic Industrial Chemicals (TICs)	Toxic Gases	cyanides, chlorine, phosgene, sulfur dioxide	BB Demo
	Toxic Liquids	BTEX, HCl, chloroform	Limited Data/TBD
	Pesticides	organophosphates/chlorines, carbamates, pyrethroids	BB Demo
	Fumigants	methyl bromide, phosphine, formaldehyde	Limited Data/TBD
	Refrigerants	HFCs, HCFCs	BB Demo
	Fire Safety Standards	PAHs	Limited Data/TBD
Explosives	Nitroaromatics	TNT, tetryl	BB Demo
	Nitramines & Nitrate Esters	RDX, PETN, HMX, EGDN, NG, Semtex, Detasheet, C-4	BB Demo
	Inorganic Nitrates	UN, AN, black powder, nitric acid	BB Demo
	Chlorate/Perchlorate	KClO ₃ , KClO ₄ , NaClO ₃ , NaClO ₄	Limited Data/TBD
	Peroxides	TATP, HMTD, MEKP	BB Demo
Drugs of Abuse	Stimulants	cocaine, amphetamine, MDMA	BB Demo
	Depressants	barbiturates, GHB	BB Demo
	Hallucinogens	LSD, psilocybin, dimethyltryptamine	BB Demo
	Opioids	heroin, codeine, morphine, fentanyl, oxycodone	BB Demo
	Cannabis	marijuana related products	BB Demo
	New Psychoactive Substances	synthetic cannabinoids/cathinones, ketamine, khat	Limited Data/TBD

BB Demo Limited Data/TBD Not Detected


Example: Chemical Warfare Agents

The vapor sensitivity of the next-generation HPMS device lies at the low parts-per billion by volume (ppbv) level with response times on the order of a few seconds. Below is an example of direct vapor detection of Sarin (GB) and thermal desorption of VX residue.

Sarin (GB) Direct Vapor Analysis




VX Thermal Desorption




Example: Explosive Residues

Negative mode MS enables the detection of several organonitrate and some inorganic explosive compounds. An integrated thermal desorber provides controlled heating of a sample collected on a swab. Trace (<500 ng) detection of pure PETN as well as detection of both PETN and RDX in a Semtex H formulation is shown below.




Example: Narcotics / Drugs of Abuse


Expanded drug detection capabilities at mass loads <50 ng currently include cocaine, heroin, fentanyl and its derivatives (see below), MDMA, amphetamine, methamphetamine, cathinone, ephedrine, pseudoephedrine, norpseudoephedrine HCl, LSD, GHB sodium salt, oxycodone and PCP. Here is one example of the specificity that can be achieved with a single compound."




fentanyl (336.5 amu)




remifentanyl (376.5 amu)



3-methylfentanyl (350.5 amu)



carfentanyl (394.5 amu)



Conclusions

The fast-switching, dual-polarity of this next-generation HPMS device, offers true trace detection capabilities in both vapor (LODs in the 10s of ppbv) and thermal desorption (LODs ~10s of ng on swab) modes. This platform is a comprehensive multi-mission solution that analyzes chemical warfare agents (CWA), the full range of explosives, high priority drugs, and relevant precursors on a single sample.

908 Devices

+1.857.254.1500 | info@908devices.com
www.908devices.com

This device is subject to export controls including those of the Export Administration Regulations of the U.S. Department of Commerce, which may restrict or require licenses for the export of product from the United States and their re-export to and from other countries. Patented technology www.908devices.com/patents
 © 2017 908 Devices