

# ZipChip

## Using the ZipChip infusion function

### OVERVIEW

The v1.3.0 release of the ZipChip software gives users the ability to use their existing ZipChip system as a static nano-spray source to infuse sample for extended ESI-MS. This function can be particularly useful for long duration MS experiments, such as top down MS of proteins, or for tuning your mass spec settings or optimizing a mass spec method.

### HOW ZIPCHIP INFUSIONS WORK

The new infusion function pushes sample into the microfluidic separation channel at the same time that the electric field is applied. So rather than forming a discreet band of sample which separates as it migrates down the channel, infusions generate a constant stream of sample for as long as the infusion is running. The infusion function can be run with any of the existing types of ZipChips but the process will be faster on High Speed (HS) ZipChips. A

screenshot of the Infusion page of ZipChip application is shown in Figure 1.

ZipChip infusions push all components of the sample matrix into the separation channel; and the infusion process does not separate analytes from salts or surfactants the way that ZipChip separations do. Therefore, it is very important that you follow the “Infusion Rules” posted on the Infuse page of the ZipChip app (Figure 1). The sample needs to be diluted into the same background electrolyte that was used to prime the chip, or else mismatches in electrical conductivity and pH could damage the chip. High concentrations of proteins can also lead to obstruction of the ESI orifice. We therefore recommend that you limit the concentration of proteins to 0.1 mg/mL or less, and promptly clear the separation channel at the end of your infusion.




Figure 1. Screen shot of the Infuse pane of the ZipChip application, version 1.3.0. The image is from a ZipChip system with an auto-sampler. Note that for Manual units, the “Rinse Sample Well” button is not available.

## EXAMPLE ZIPCHIP INFUSION: NATIVE NIST MAB ON AN HRN CHIP

The infusion process is particularly useful for optimizing mass spec settings for native protein analysis. Our application scientists run this specific infusion process whenever they attempt native antibody work on a new mass spec for the first time. Because it uses the same BGE and chip type as our ZipChip native antibody separation methods, users can switch from infusion to separation without having to re-prime a chip.

### Preparing the chip and the sample:

Begin by priming an HRN chip with Native Antibody BGE. Dilute 10  $\mu$ L of the NIST Monoclonal Antibody Reference Material with 90  $\mu$ L of LC-MS grade water. This creates a 1 mg/mL sample which can later be used for ZipChip separations. **Do not attempt to use this sample for infusions** – it is too concentrated and its sample matrix does not match the BGE in the chip. Dilute 10  $\mu$ L of the 1 mg/mL sample with 90  $\mu$ L of Native Antibody BGE. This creates a 0.1 mg/mL sample which is appropriate for ZipChip infusion. The sample must be manually loaded into the chip. Slide the chip out of the chip manifold and manually remove any liquid that is present in the sample well (see Figure 2) of the chip. Pipet 20 to 40  $\mu$ L of the 0.1 mg/mL sample into the sample well, then lock the chip back into the chip manifold.




Figure 2. ZipChip schematic showing location of the Sample well.

### Acquiring data during the infusion process:

To demonstrate all of the features of the ZipChip infusion, this example includes acquisition of MS data for the entire process. The Acquisition settings shown in Figure 3 were used. Note that a Method file was not used so that the settings could be adjusted in real time. Also note that the Acquisition time was set to “continuously”, and

“On start” was set to “don’t wait”. The ZipChip infusion process does not trigger the mass spec to start data collection the way that a separation method does, so do not set the MS to wait for a contact closure signal. The continuous acquisition means that the data acquisition must be manually stopped when finished.




Figure 3. MS Tune Acquisition settings.

### Running the Infusion:

On the Infuse page of the ZipChip app (see Figure 1), set the *Infusion Field Strength* to 500 V/cm, then press the *Start Infusion* button. The electrospray process begins immediately. For this example, the Start button on the MS Tune Acquisition was pressed immediately after starting the infusion. Figure 4 shows the base peak ion signal recorded during the infusion process. Note that it takes about 5 minutes before the antibody signal is detected by the mass spectrometer. This is the time it takes for the antibody molecules to migrate down the separation channel. Once the antibody reaches the ESI orifice, a constant mass spec signal is detected for as long as the infusion process is run. Figure 5 shows a mass spectrum from this example infusion of the NIST mAb. For this example, the Stop button on the Infuse page of the ZipChip app was pressed 15 minutes after the infusion was started and about 10 minutes after the antibody was first detected. If necessary the sample can be infused for longer times, but it is recommended to refresh the BGE and replace the sample every hour. As seen in Figure 4 the MS signal immediately dropped to zero when the Stop button was pressed since the electrospray process stopped. At that point the entire separation channel was still full of sample, so the *Clear Channel* button was immediately pressed to begin emptying the channel. It took an additional 6 -7 minutes for the residual antibody to fully clear out of the separation channel. At this point the microfluidic channels of the chip are clean. To proceed with ZipChip separations or another infusion experiment, simply rinse the sample well and refresh the BGE.




Figure 4. Mass spec base peak intensity observed during the example infusion of the NIST mAb on an HRN chip.




Figure 5. Mass Spectrum from NIST mAb infusion. The inset is a zoomed view of the most intense charge state.

## SUMMARY

The ZipChip system can now be used as a static nano-spray source for infusing samples. The new Infuse pane in the ZipChip software provides a quick and easy way to infuse samples with the same hardware and consumables used for ZipChip separations. ZipChip infusions can

be used for a wide variety of applications such as top-down MS/MS fragmentation, ultra-high-resolution MS characterization, or optimization of MS settings for new assays. This feature combined with ZipChip separations makes the system a powerful and versatile tool well suited for in-depth characterization of complex samples.



**908devices**  
www.908devices.com

ZIPCHIP IS FOR RESEARCH USE ONLY

Zip Chip is subject to export controls including those of the Export Administration Regulations of the U.S. Department of Commerce, which may restrict or require licenses for the export of product from the United States and their re-export to and from other countries. Patented technology www.908devices.com/patents © 2019 908 Devices