Microchip CE-MS Analysis of Sialylated N-glycan linkage Isomers Mengxia Cheng¹, Haojie Lu^{1, *}, Liang Wang², Kate Yu²



Overview

- Linkage isomers of sialylated N-glycans are of great significance for diagnosing and monitoring diseases.
- Qualitative and quantitative analysis of sialylated Nglycan linkage isomers remains challenging due to their low abundance and limited isomeric separation techniques.
- A novel strategy based on microchip capillary electrophoresis-mass spectrometry (Microchip CE-MS) was developed.
- Fast and specific analysis of **α-2,3- and α-2,6-linked** sialylated N-glycan linkage isomers can be achieved.

Methods

- A novel strategy based on CE-MS workflow was utilized for this work combining fast CE separation with high-resolution mass spec identification.
- Front-end CE and NanoESI was performed by ZipChip[™] interface from 908Devices Inc.
- Only positively charged components in the sample can be analyzed based unique microchip CE principle.
- CE separation was implemented in a liquid BGE solution composed of the mixture of methanol and water (ZipChipTM Metabolite assay kit) through a 22 cm length CE separation channel (ZipChip[™] HR chip).
- MS analyses were performed by interfacing the device with a QExactive HF MS from Thermo Fisher.



- N-glycans from glycopeptides or glycoproteins were prepared with standard enzymatic digestion procedure (PNGase F in ABC buffer).
- N-glycans from serum samples were prepared through quantification (Protein Assay Kit), reduction (DTT), alkylation (IAA), followed by enzymatic digestion with PNGase F.
- Released Glycan were further derivatized on the sialic acid site, the different reagents were used depend on the application. The solution was then purified by cotton-wool-based HILIC-SPE.

^{1.}Department of Chemistry and Institutes of Biomedical Sciences, Fudan university, Shanghai, 200032, China; ^{2.}908Devices Inc. Boston, MA, 02210, USA



