# 

ON-LINE MONITORING AND REAL-TIME CONTROL OF GLUCOSE AND LACTATE



# MEET MAVEN

# OPTIMIZED, AUTOMIZED, DON'T-EVEN-SWEAT-IT GLUCOSE AND LACTATE CONTROL.

Monitoring and control of glucose and lactate concentrations during cell culture and fermentation are critical to achieve optimal process efficiency and productivity for effortless realization of PAT strategies. Real time insights into cellular metabolism are invaluable for media optimization and development of feeding and control strategies.





Detect and react to process changes in real time with on-line monitoring of glucose and lactate as often as every 2 minutes.



Develop advanced process control strategies and identify opportunities for process improvement with up to 720 glucose-lactate data points per day.



Preserve precious bioreactor material and reduce process contamination risk with sample-free analysis based on the proprietary selective flow diffusion method.



Reliable and precise biosensor-based measurement technology ensures that optimal glucose concentration can be maintained, even at very low levels.



Reusable and single-use probes and flow cells make it easy to integrate MAVEN with most types and sizes of upstream systems in both the PD lab and GMP manufacturing.




Control the feeding pump directly with the integrated PID controller or easily connect with existing process automation systems via OPC UA, Analog or Serial outputs.

### DESIGNED WITH FEATURES TO MAKE YOUR DAY

KNOW EXACTLY WHAT'S HAPPENING IN YOUR BIOREACTOR, MAINTAIN CONTROL OF YOUR PROCESS, AND KEEP YOUR CELLS HAPPY.

Unlike in manual or automated sampling, measurements are performed without any loss of bioreactor volume and without increased risk of process contamination.


Glucose and lactate molecules from the medium diffuse through the semi-permeable membrane into the transport buffer solution circulated between the probe and measuring cell.





#### HIGHLY PRECISE AND ROBUST BIOSENSORS

The biosensor's enzymatic recognition and selectivity for glucose and lactate generates an electrical signal proportional to concentration of targeted analytes.



- Sensor type:Enzyme-based electrochemical biosensor
- Glucose measurement range:
   From 0.01 to 40 g/L concentration
- Lactate measurement range:From 0.05 to 10 g/L concentration
- Replaceable biosensor service life:
   5,000 measurements or 30 days

#### BE A REBEL AND REAP THE REWARDS


Combine REBEL's at-line analysis of amino acids and vitamins with MAVEN's on-line monitoring and control of glucose and lactate to get unambiguous answers where and when you need them.





## HOW WILL YOU USE MAVEN?


MONITOR GLUCOSE AND LACTATE IN YOUR BIOREACTOR IN REAL-TIME. REFINE FEEDING STRATEGIES ON THE FLY.



### MAVEN LOVES TO COLLABORATE

LIKE A GOOD LAB PARTNER, MAVEN EASILY INTEGRATES WITH YOUR EXISTING OR NEW UPSTREAM PROCESS DEVELOPMENT OR GMP MANUFACTURING SYSTEMS.

✓ Compatible with Most Bioreactors and Fermenters
✓ Automation Included
✓ Advanced Connectivity



### SUDDENLY, YOU'RE THE BOSS OF YOUR BIOREACTOR


Tight control of glucose-lactate concentrations enables development and implementation of dynamic feeding strategies to optimize cellular metabolisms leading to higher productivity and improved quality of biotherapeutics.

Bolus Control of Feeding


> Dynamic Control of Feeding



Maintain lower range of glucose concentration and minimize nutrient fluctuations



Optimize lactate metabolism and prevent growth-inhibiting lactate accumulation



Improve cell growth, increase maximum viable cell density, and total product yield



Improve product quality such as glycosylation profile and reduce unwanted glycation

**≫908** devices

### MAVEN

#### MAVEN@908DEVICES.COM

MAVEN is subject to export controls including those of the Export Administration Regulations of the U.S. Department of Commerce, which may restrict or require licenses for the export of product from the United States and their re-export to and from other countries. Patented technology 908devices.com/patents/© 2022 908 Devices. LS-MN BR-0001 01-2023