

2E8

1E7

Online Peptide Fractionation

Can we modify the ZipChip-SPE-CE-MS method to elute the peptides off of the SPE bed in multiple fractions?

Particles = 5 μ m porous silica, C18 **BGE = Acetonitrile/Water/Formic acid (ZipChip Peptides BGE)** Weak mobile phase (WMP) = 1% acetonitrile, 0.1 M ammonium acetate

Sample: • 20 ng/µL HeLa Tryptic digest

Method:

- ZipChip SPE sample loaded as normal (240 seconds)
- Run manually stopped after the Wash step
- Sequence of Standard ZipChip runs performed with increasing % acetonitrile elution solvents in place of samples

Highlights:

- Successfully eluted different peptides in different fractions
- Total IDs from the fraction runs were ~2x more than running a single ZipChip-SPE run of the same sample

How to improve:

- Adjust acetonitrile content of fractions to more evenly distribute the peptides
- Include a WMP wash before running each fraction • This will yield sharper peaks for fractions 2+

Size Exclusion Desalting for Charge Variant Analysis of mAbs

Instead of solid phase extraction, could we pack the bed with SEC particles and use it for online desalting of Biotherapeutic Proteins

Particles = 5 μ m porous silica, Diol **BGE = Acetonitrile/Water/Formic acid (ZipChip Peptides BGE)**

Sample:

• NIST mAb diluted from stock formulation to 0.25 mg/mL with 1x PBS

Method:

- Standard ZipChip method with injection volume setting increased to compensate for resistance of packed bed
- Pressure assist on to back flush salt out of SPE bed during each run

Highlights:

- Effective online desalting of mAb
- Well resolved charge variant peaks from a salty sample

How to improve:

- Use better particles
- Tighter size distribution and smaller pores
- Make some adjustments to the ZipChip method to achieve more reproducible performance of this application



Exploring the Possibilities for Microchip SPE-CE-MS





If we swap the C18 particles for more protein-friendly particles, can we use the ZipChip-SPE-CE-MS method for analysis of intact proteins?

Particles = 5 μ m porous polymeric particles **BGE = Acetonitrile/Water/Formic acid (ZipChip Peptides BGE)** WMP = 1% acetonitrile, 0.1 M ammonium acetate

Method:

• Standard ZipChip SPE-CE-MS method

Highlights:

- The proteins in the Pierce Intact Protein mix behaved well with this method
- Concentration factor >100x compared to regular ZipChip without tITP focusing
- Concentration factor >20x compared to regular ZipChip with tITP focusing

How to improve:

- Automate a more thorough rinse of the bed between runs to minimize carryover
- Adjust the solvent compositions to optimize performance for intact proteins

Step 2: Sample in 0.1 M AmAc, up to 3% MeCN, various carriers/additives

Online SPE Desalting of mAbs for Bioreactor Monitoring

Could we use a version of the Top-down method for simple desalting of mAbs pulled directly from a bioreactor?

Particles = 5 μ m porous polymeric particles BGE = Acetonitrile/Water/Formic acid (ZipChip Peptides BGE) WMP = 1% acetonitrile, 0.1 M ammonium acetate

Samples:

- NIST mAb diluted from stock formulation with CDCHO
- Each sample was further diluted 10x with ZipChip SPE Diluent
- Concentrations listed were before 10x dilution with diluent

Method:

Standard ZipChip SPE-CE-MS method

Highlights:

- Effective desalting and concentration of mAb
- Clean spectra with accurate deconvolution at 10 μ g/mL in CDCHO

How to improve:

- Automate a more thorough rinse of the bed between runs • Adjust the solvent compositions to optimize performance for intact
- Explore a combined method that measures unretained small molecules injected through the bed on the same device

Top Down Proteomics





The technologies discussed in this poster are the subject of one or more granted/pending patents. www.908devices.com/patents/