

Overview

Rapid charge variant separation of mAbs and ADCs by microchip capillary electrophoresis coupled with trapped ion mobility mass spectromery.

Conformational changes and alterations of average drug-to-antibody ratio were observed among charge variants along with changes in the collisional cross section between mAbs and ADCs.

Introduction

therapeutics have become one of the largest drug classes recently, fueled by the therapeutic efficacy of monomAbs) and antibody-drug conjugates (ADCs). However, mAb and ADC drug products are heterogeneous tible to formation of charge variants that may increase the rate of proteolytic degradation, alter conformation and or likelihood of aggregate formation¹. For the first time, we present a native charge variant analysis by Microchip capillary electrophoresis (CE)² coupled with trapped ion mobility spectrometry³ and mass spectrometry (CE-TIMS-MS) for four antibodies: NISTmAb, a native cysteine-linked ADC and its parent mAb, and lysine-linked ADC (Kadcyla). The developed method measures gas phase collisional cross section (TIMSCCS_{N2}) for conformational state assessment of individual charge variants in only 10 minutes time. Our developed method enables rapid TIMSCCS_{N2} measurement for individual mAb and ADC species to assess potential conformation changes. This method revealed significant changes in charge variant previous spectroscopio



investigations of the impact of mAb deamidation on conformation^{4,5}. This method also partially separates DAR species of the lysine-linked ADC for MS analysis of the extremely complex drug products, a feat not possible by previously developed direct infusion TIMS-MS³. The native CE-TIMS-MS platform improves MS signal-to-noise relative to direct infusion, eliminates the need for sample clean up prior to injection, and uses <2 ng sample per analysis. This platform determines quality attributes such as intact mass, TIMSCCS_{N2}, and average DAR for individual charge variant species in only 10 minutes.

Workflow Acidic 2 Time [min] Native Microchip **CE-TIMS-MS**

ed charge variants for individual assesment of conformation by TIMS to enable charge variant specific TIMS CCS_{N2} information to be obtained.

Methods

Charge variants of a model cysteine-linked ADC, model lysine-linked ADC, and their respective parent mAbs were separated by a ZipChip Microchip CE Device (908 Devices) coupled with a Bruker timsTOF Pro for TIMSCCS_{N2} and intact mass measurement. Using a High Resolution Native chip and Charge Variant TOF background electrolyte, charge variants of the ADCs and parent mAbs were separated with a field strength of 500 V/cm and a 10-minute total runtime. 1 nL of sample was pressure injected, with ADCs at a concentration of 2 mg/mL and the mAbs at 1 mg/mL. The use of TIMS allowed calculation of TIMS CCS

for all mAbs and ADCs at the 26⁺ charge state using the Mason-Schamps equation (Equation 1)⁶. Where μ is the reduced mass of the ion, k_b is Boltzmann's constant, T is the drift region temperature, z is the ionic charge, e is the charge of an electron, N

the buffer gas density, and K_o is the reduced mobility. Data analysis was performed using Compass DataAnalysis 5.3 and MASH Native v1.1. The average DAR for individual regions was determined by Maximum Entropy deconvolution with a resolution of 10,000 followed by Sum Peak peak picking to calculate the area under the curve for all observed DAR species (Equation 2). Where the area under the curve for all observed DAR species weighted by drug load is divided by the total unweighted area under the curve for all DAR species. For the model cysteine-linked ADC, odd numbered drug loading was not observed and the contribution of odd numbered DAR

Equation 1 TIMS CCC

species was not included in the av- **Equation 2** Average DAR =erage DAR calculation.

[(DAR1 * 1) + (DAR2 * 2) + (DAR3 * 3) ...[(DAR0) + (DAR1) + (DAR2) + (DAR3) ...

 $LO_{N2} =$

 $| 2\pi$

 N_0K_0

Conformation of Native Antibody-Drug Conjugate Charge Variants Revealed by Microchip Capillary Electrophoresis Coupled with Trapped Ion Mobility

Eli J. Larson¹, Adi M. Kulkarni², Jake A. Melby¹, Matthew Fischer¹, Zhan Gao¹, Kevin M. Buck¹, Melissa R. Pergande¹, Yanlong Zhu¹, Guillaume Tremintin³, Kate Yu², Ying Ge¹ ¹University of Wisconsin-Madison, ²908 Devices Inc., ³Bruker Daltonics





Figure 1. Native microchip CE-TIMS-MS analysis of model therapeutic mAb in histidine buffer at pH 7 with 6% sucrose. The total ion electropherogam (A) shows three mAb peaks: the fully lysine clipped mAb (Main, 0K) and two deamidated acidic charge variants (Acidic 1 and Acidic 2). A heatmap of mobility and *m*/*z* for the main peak shows limited dependence of mobility on charge state (B). The exctracted ion mobilograms for each peak (C) shows a shift in the mobility of charge variant peaks which is also reflected in the calculated TIMSCCS_{N2} for each region (D), with the Acidic 1 variant showing a significantly larger TIMSCCS_{N2} than the main peak (p=0.002). The results support previous spectroscopic studies of the impact of deamidation on mAb conformation^{4,5}.





Acidic 2 All Peaks Acidic 1

Figure 4. Impacts of charge variant species on measured TIMSCCS_{N2} and average DAR value. Comparison of the measured TIMSCCS_{N2} for the model therapeutic mAb and ADC revealed a significant difference (p=0.009) in the ™ ^sCCS_{No} of the first acidic variant (A). The calculated aver-3 age DAR values for charge variants of the ADC (B) showed significant changes in the average DAR for all comparison points (Main-Acidic 1: p=0.003; Main-Acidic \checkmark 2: p=0.001; Acidic 1-Acidic 2: p=0.006). The average DAR calculated across all charge variant peaks aligns with previously determined values³. Plotting the TIMSCCS_{N2} and average DAR for each charge variant of the model cysteine-linked ADC (C) underscores the opposing trends observed: deamidation variants show an increase in CCS, but decrease in average DAR.



tropherogam (A) shows that charge variants, whether caused by modifications or drug conjugation events, are not fully resolved due to the extreme complexity of the sample. A mobility vs m/z heatmap from 8.6-8.8 minutes migration time (B) show the heterogenity of overlapping charge states and overall low signal for unique DAR species due to drug product heterogenity. The charge deconvoluted MS1 spectra of ten 0.2 minute migration regions (C) demonstrates that as migration time increases, high DAR species (DAR 4-7) are predominantly observed. Considering the average DAR for each migration region reveals a very strong linear correlation between migration time and average DAR indicating partial separation by level of drug conjugation, though DAR species are unresolved (D). The TIMSCCS_{N2} calculated for each migration region exhibits larger standard deviations than those observed for the mAbs and cystine-linked ADC, potentially due to coelution of multiple charge variant species (E).

Conclusion

The combination of microchip capillary electrophoresis and trapped ion mobility enable a rapid, informative method for native monoclonal antibody and antibody-drug conjugate charge variant analysis. Our method allowed detection of significant TIMSCCS_{N2} changes between the unmodifed main peak and deamidated peaks, and can be used to identify charge variants that pose a threat to overall drug stability. Leveraging the powerful combination of a rapid, native liquid-phase separation with TIMS, our approach provides a new tool for pre- and post-clinical quality control which can serve as a foundation to address the ever more complex array of antibody-based therapeutics currently under development.

References

- 5) Lu, X. et al. Sci. Rep. 2020, 10, 383.
- Gabelica, V. et al. Mass Spectrom. Rev. 2019, 38, 291-320
- /) Upton, R. et al. Chem. Sci. 2019. 10, 2811-2820.
- 2) Wu, Z. and Wang, H. et al. J. Pharm. Biomed. Anal. 2023, 223, 11514

1) Khawli, L. A. et al. mAbs. 2010, 6, 613-624

- 3) Larson, E. J. and Roberts, D. S. et al. Anal. Chem. 2021. 93, 10013-10021
- 4) Lu, X. et al. mAbs. 2018. 11, 45-57.



Acknowledgements

The authors acknowledge the support of NIH grants R01HL109810-05A1 R01GM125085, R01HL109810, and R01HL096971.





