

In-depth Characterization of Adeno-Associated Viruses (AAVs) using Microchip CE-MS

Josh Smith¹, Sara Carillo¹, Adi M Kulkarni², Erin Redman³, Kate Yu², Jonathan Bones^{1,4}

¹Characterization and Comparability Laboratory, NIBRT, Dublin, Ireland; ²908 Devices, Inc., Boston, MA; ³908 Devices, Inc., Carrboro, NC; ⁴School of Chemical and Bioprocess Engineering, University College Dublin, Belfield, Dublin, Ireland



Introduction

- Adeno-associated viruses (AAVs) are promising vectors for gene therapy products due to their low immunogenicity.
- These icosahedral, 60-mer capsids are composed of the three viral proteins (VPs) VP1, VP2, and VP3, in an approximate 1:1:10 ratio, respectively.
- Characterization of AAV viral capsid proteins (VPs) and their post translational modifications (PTMs) is critical to ensure product quality and safety.
- Commonly utilized liquid chromatography–mass spectrometry (LC-MS) characterization methods are often time consuming, with analysis alone often taking 30 minutes to >2 hours, depending on the desired characterization application.
- AAV sample is often limited, so minimizing the sample required for analysis is crucial.
- Microchip capillary electrophoresis coupled with MS (CE-MS) addresses such concerns by providing rapid characterization capabilities with minimal sample. It has previously proven capable of performing peptide mapping on AAVs in under 10 minutes [1] and intact analysis of capsid VPs in under 5 minutes [2].
- Here we demonstrate the utilization of microchip CE-MS for VP characterization and peptide mapping for rapid in-depth AAV characterization. Performed across a spectrum of AAV serotypes, we illustrate how this platform quickly provides crucial information related to the characteristics of AAVs including the presence of VP proteoforms and PTMs.

Viral Capsid Protein (VP) Analysis

Methodology

VP Separation

1 in 5 dilution with Peptide BGE

Thermomixer

15 min, 37°C, 500 rpm

ZipChip-MS Analysis

OE240

Data Processing


BioPharma Finder™ 5.1

Acquisition Settings

Field Strength: 500 V/cm	Injection Volume: 5.5nL
Application Mode: Intact protein, Pressure Mode: Low Pressure, Expected Peak Width: 10	Default Charge State: 35
MS1	15k Res, m/z 740-2000, RF Lens 125%, Normalized AGC 50%, max IT 200 ms, Collision Energy: 35V

Limit of Detection (LOD) Analysis



Figure 2: LOD testing of AAVs for microchip CE-MS using Empty AAV8. All Identified VP proteoforms were detected at concentrations ≥8.80E6 viral particles, with at least one proteoform of each VP detected at concentrations as low as 2.64E6 viral particles



Viral Capsid Protein (VP) Analysis Continued

Proteoform Identification

Figure 3: Illustration of reproducibility of VP analysis performed using microchip CE-MS



Peptide Mapping

Methodology

Figure 6: Workflow for AAV Peptide mapping using microchip CE-MS

Peptide Mapping Acquisition Settings


ZipChip	Field Strength: 500 V/cm, Injection Volume: 8.0nL
Global Settings	Application Mode: Peptides, Pressure Mode: Standard, Expected Peak Width: 15, Default Charge State: 2
MS1	Polarity: Positive, 30K Res, m/z 200-1800, RF Lens: 50%, Normalized AGC: 100%, Microscans: 1, Max IT: 200 ms
MS2	Isolation window (m/z): 2, NCE: 28, 15K Res, Ref Lens: 70%, AGC 100%, Microscans: 1 Max IT: 50 ms
DDA	Cycle Time, Scan Event: 1, Dynamic Exclusion: 1s Mass Tolerance: 10ppm, Charge State: 2-8, Intensity Threshold: 5.0e4

Samples: Same as for VP Separation

Sample Preparation: 10 µg of AAV (determined using NanoOrange protein assay) was digested SMART Digest pepsin magnetic beads on a Kingfisher™ DuoPrime. Desalting is performed using a C18 spin column and then reduced to dryness in a speed-vac.

Sample Analysis: Separation is performed using a ZipChip™ Device coupled to an Orbitrap Exploris™ 240 MS. Sample is reconstituted in 10uL of Peptide BGE and loaded directly into the well of an HR chip primed with Peptide BGE for analysis. Each run is 20 minutes, and 2 injections were performed for each sample.

Software: Same as for VP Separation



Conclusions & Future Work

- Rapid microchip CE-MS can be utilized for the in-depth characterization of AAVs including viral capsid protein analysis and peptide mapping.
- VP separation is performed in as little as 5 minutes while detecting multiple proteoforms of each viral capsid protein.
- Peptide mapping provides high sequence coverage of all VP proteins in under 20 minutes and can characterize PTMs.
- Application demonstrated with serotypes AAV6, AAV8, and AAV9 to show the broad applicability of this platform across AAV serotypes.
- With characterization performed at greatly increased speeds compared to standard LC-MS analysis, the microchip CE-MS ZipChip platform can provide an attractive alternative for laboratories looking for high throughput AAV characterization.
- Future work:** Investigate VP ratios; Improve detection of PTMs on N-terminus of VP1; Apply microchip CE-MS to empty/full AAV intact capsid analysis