

Characterization of Protein/PTM Changes by Chip Capillary Electrophoresis Mass Spectrometry

Deborah Chang; Deborah R Leon; Christian F Heckendorf; Joseph Zaia; Catherine E. Costello; Mark E. McComb


Boston University School of Medicine, Boston MA

BOSTON
UNIVERSITY

Overview

- Characterization of post-translational modifications (PTMs), specifically complex protein glycosylation, remains a challenge for high-throughput mass spectrometry applications.
- Standard C18 nano-LCMS approaches are sufficient for certain PTM compound classes but are unable to adequately resolve individual glycoforms.
- Capillary electrophoresis mass spectrometry (CE-MS) offers high-resolution separation capabilities which are unique and potentially offer advantages for glycoproteomics analyses.
- We successfully conducted a series of experiments in collaboration with 908 Devices (Boston, MA) using chip-based CE-MS and demonstrated the ability to obtain high-resolution separation on digests of proteins and glycoproteins.(1)
- Here we present results of this continued work and demonstrate improvements in the entire analytical pipeline including sample preparation, separation parameters, MS and MS/MS data acquisition and bioinformatics approaches for data interpretation.

Applications of CEMS




Approach

- CE was performed using both prototype and commercially available second generation ZipChip™ CE interfaces obtained through collaboration with 908 Devices.
- CE-MS analyses were performed by interfacing the device with a QExactive MS or a QExactive HF MS (Thermo-Fisher).
- Samples included in-solution protein digests with emphasis on protein glycosylation.
- Peptides/proteins were prepared in up to 50% organic solvent (ACN, MeOH or isopropanol) prior to CEMS; picomole to low femtomole sample loads were analyzed.
- MS/MS was performed using DDA and AIF modes which were optimized for high-throughput bottom-up acquisition. A wide range of MS, MS/MS and AIF parameters were explored.
- Feature identification was enabled by analyzing the data using Proteome Discoverer (Thermo-Fisher), Mascot (Matrix Science), Scaffold (Proteome Software), PEAKS Studio (Bioinformatics Solutions) and Byonic (Protein Metrics).
- Detailed assignments of glycopeptides were achieved using GlycReSoft(2) and BUPID-PALM (Boston University Protein Identifier-Parallel Acquisition Library Matcher) software developed in our laboratory.
- Top down MS and MS/MS analysis was performed using BUPID software developed in our laboratory at BUSM.

Results


- While we have demonstrated that chip CE-MS allows for unique separation of glycopeptides (1), we have observed that we sometimes obtain broad peaks due to the presence of sialic acid termini on N- or O-linked glycans.
- We also observed that, in some cases, due to fast separation and narrow peaks, the quality of MS/MS data was relatively poor or precursors were missed.
- We further explored different electrolytes to improve CE separation, evaluated different MS and MS/MS data acquisition parameters to accommodate narrow peak widths, and expanded upon our initial set of CE and MS parameters.
- Chip-based CE-MS of proteolytic digests of A1AG1 yielded 7-minute to 15-minute run times, depending on conditions and obtained data comparable with >60-minute LC-MS experiments.

Comparison of CE-MS with LC-MS: A1AG1 Peptides



- Setting the m/z range 600-3000 resulted in an increase in S/N compared with m/z 370-1800 and did not result in the loss of glycopeptide MS/MS data.
- MS2 threshold was set at 2e5 to 5e5 and yielded more confident assignments of higher quality MS2 data.
- Reducing the duty cycle to Top 4 allowed for more cycle time, improved MS1 TIE (SIE) and accommodated narrow CE peaks.
- Glycopeptide MS/MS required 100msec fill/transient at 17,500 resolution. We observed that 2 scans yielded improved glycopeptide MS/MS data.

CE-MS: A1AG1 Peptides/Asialo-Glycopeptides




- Use of different enzymes allowed for improved digestion to peptides, and more control on peak capacity.
- Desialylation improved peak shapes and CE run times.
- Arg-C digests yielded larger glycopeptides which behaved well under the CE conditions of separation.


Results

- Glycopeptides were consistently well resolved from peptide peaks by CE-MS.
- The individual glycopeptides migrated in terms of increasing peptide mass and oligosaccharide composition.
- This improvement in resolution compared with LCMS would potentially allow for improved differential glycan structure quantification.

CE-MS of A1AG1 Glycopeptides: Increasing Migration Times




Migration Times: CE-MS Chymotrypsin Peptides FYFTPNKTEDTIF



- Optimal glycopeptide MS/MS data requires longer duty cycles in comparison to peptide MS/MS: 100 msec and 2 scans. This resulted in loss of identifications for lower S/N glycopeptides.
- Use of targeted exclusion and inclusion lists via a 2-run strategy improved sequence coverage and assignment of glycopeptides.
- Implementation of wide precursor windows, AIF and alternating high/low collision energy also improved our ability to obtain high quality sequence information for glycopeptides.
- In order to interpret data obtained using large precursor window isolation, AIF and high/low fragmentation modes, we utilized a spectral library search approach.


A1AG1: CE-MS AIF MS/MS: Glycopeptide ID via BUPID-PALM



Results

- Lastly we explored the use of Top Down MS and MS/MS for glycosylation characterization.
- Intact RNase-B was used as a model glycoprotein due to its low molecular weight and single glycosylation site.
- Different CE voltages and electrolytes were used.
- MS and MS/MS parameters were varied in a similar means as the bottom up experiments.
- We also implemented targeted MS/MS windows and AIF modes for acquisition.
- We were able to separate RNase-B from RNase-A.
- Resolution of the different high-mannose forms of RNase-B was achieved.
- Both targeted MS/MS and AIF modes yielded sufficient information to characterize RNase-B and the individual glycoforms.

Top Down CE MS/MS of RNase-B



Summary

- These strategies further improved the quality of the data and allowed for a greater number of confident glycopeptide assignments.
- In conclusion we demonstrate improved approaches for CE-MS for characterization of complex protein glycosylation.
- See poster MP 359 BUPID-PALM: Glycopeptide Identification by All-Ion Fragmentation(AIF) Ion Mobility MS/MS; Christian Heckendorf; Joshua A. Klein; James A. Hill; Catherine E. Costello; Mark E. McComb; Boston University School of Medicine, Boston, MA

Acknowledgements

NIH grant P41 GM104603
NHLBI Contract HHSN268201000031C
We thank 908 Devices for access to the ZipChip CEMS Interface and J. Michael Ramsey for helpful advice.

We thank all contributing members of
the Center for Biomedical Mass Spectrometry

- (1) Khatri K, et al. Anal. Chem. 2017 PMC554952.
- (2) Maxwell, E., et al. PLoS ONE 2012 PMC3458864.