

Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law Issue date: 8/2/2023 Revision date: 10/26/2023 Version: 1.1

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Mixture

Product name : ZipChip Peptides Diluent

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category : Professional use

Use of the substance/mixture : For research and development use only.

Function or use category : Laboratory chemicals

1.2.2. Uses advised against

Restrictions on use : Not for use in diagnostic procedures.

1.3. Details of the supplier of the safety data sheet

Manufacturer

908 Devices 645 Summer St 02210 Boston, MA USA

T 1 (857) 254 - 1500 <u>908devices.com</u>

1.4. Emergency telephone number

Emergency number : 1 (844) 908 - 4357

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

 Flam. Liq. 2
 H225

 Acute Tox. 4 (Oral)
 H302

 Eye Irrit. 2
 H319

Full text of hazard classes, H- and EUH-statements: see section 16

Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP) :





GHS02

GHS07

Signal word (CLP) : Danger

Contains : Acetonitrile; Formic acid

Hazard statements (CLP) : H225 - Highly flammable liquid and vapour.

H302 - Harmful if swallowed.

H319 - Causes serious eye irritation.

Precautionary statements (CLP) : P210 - Keep away from heat, hot surfaces, sparks, open flames and other ignition sources.

No smoking.

P240 - Ground and bond container and receiving equipment.

P280 - Wear protective gloves/protective clothing/eye protection/face protection/hearing

Safety Data Sheet

Unknown acute toxicity (CLP) - SDS

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

protection

P337+P313 - If eye irritation persists: Get medical advice/attention.

P403+P235 - Store in a well-ventilated place. Keep cool.

P501 - Dispose of contents/container to hazardous or special waste collection point, in

accordance with local, regional, national and/or international regulation.

1% of the mixture consists of ingredient(s) of unknown acute toxicity (Oral)

1% of the mixture consists of ingredient(s) of unknown acute toxicity (Dermal)
1% of the mixture consists of ingredient(s) of unknown acute toxicity (Inhalation (Vapours))

Unknown hazards to the aquatic environment (CLP) : Contains 1 % of components with unknown hazards to the aquatic environment

2.3. Other hazards

Contains no PBT/vPvB substances ≥ 0.1% assessed in accordance with REACH Annex XIII

SECTION 3: Composition/information on ingredients

3.1. Substances

Not applicable

3.2. Mixtures

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
Acetonitrile substance with national workplace exposure limit(s) (GB); substance with a Community workplace exposure limit	CAS-No.: 75-05-8 EC-No.: 200-835-2 EC Index-No.: 608-001-00-3	< 50	Flam. Liq. 2, H225 Acute Tox. 4 (Oral), H302 (ATE=500 mg/kg bodyweight) Acute Tox. 4 (Dermal), H312 (ATE=1100 mg/kg bodyweight) Acute Tox. 4 (Inhalation), H332 (ATE=11 mg/l/4h) Acute Tox. 4 (Inhalation:vapour), H332 (ATE=11 mg/l/4h) Eye Irrit. 2, H319
Formic acid substance with national workplace exposure limit(s) (GB); substance with a Community workplace exposure limit (Note B)	CAS-No.: 64-18-6 EC-No.: 200-579-1 EC Index-No.: 607-001-00-0	< 2	Flam. Liq. 3, H226 Acute Tox. 4 (Oral), H302 (ATE=1100 mg/kg bodyweight) Acute Tox. 3 (Inhalation:vapour), H331 (ATE=7.85 mg/l/4h) Skin Corr. 1A, H314 Eye Dam. 1, H318

Specific concentration limits:		
Name	Product identifier	Specific concentration limits
Formic acid	EC-No.: 200-579-1 EC Index-No.: 607-001-00-0	(2 ≤ C < 10) Skin Irrit. 2, H315 (2 ≤ C < 10) Eye Irrit. 2, H319 (10 ≤ C < 90) Skin Corr. 1B, H314 (90 ≤ C < 100) Skin Corr. 1A, H314

Note B - Some substances (acids, bases, etc.) are placed on the market in aqueous solutions at various concentrations and, therefore, these solutions require different classification and labelling since the hazards vary at different concentrations. In Part 3 entries with Note B have a general designation of the following type: '... %'. In this case the supplier must state the percentage concentration of the solution on the label. Unless otherwise stated, it is assumed that the percentage concentration is calculated on a weight/weight basis.

Full text of H- and EUH-statements: see section 16

Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures after inhalation : If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable

for breathing. Get medical advice/attention if you feel unwell.

First-aid measures after skin contact : IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with

water/shower. Wash clothing before re-using. Get medical attention if irritation develops and

persists.

First-aid measures after eye contact : IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if

present and easy to do. Continue rinsing. If eye irritation persists: Get medical

. advice/attention

First-aid measures after ingestion : IF SWALLOWED: rinse mouth. Do NOT induce vomiting. Call a POISON CENTER/doctor if

you feel unwell. Never give anything by mouth to an unconscious person.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/effects after inhalation : May cause irritation to the respiratory tract.

Symptoms/effects after skin contact : May cause skin irritation. Repeated exposure may cause skin dryness or cracking.

: Causes serious eye irritation. Symptoms may include discomfort or pain, excess blinking

and tear production, with marked redness and swelling of the conjunctiva.

Symptoms/effects after ingestion : Harmful if swallowed. May cause gastrointestinal irritation, nausea, vomiting and diarrhea.

4.3. Indication of any immediate medical attention and special treatment needed

Symptoms may be delayed. In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

SECTION 5: Firefighting measures

5.1. Extinguishing media

Symptoms/effects after eye contact

Suitable extinguishing media : Dry chemical powder. Alcohol resistant foam. Carbon dioxide (CO2). Unsuitable extinguishing media : Do not use a solid water stream as it may scatter and spread fire.

5.2. Special hazards arising from the substance or mixture

Fire hazard : Highly flammable liquid and vapour. Products of combustion may include, and are not

limited to: oxides of carbon. Nitrogen oxides. Hydrogen cyanide. Ammonia. Acetic acid.

Toxic vapours. irritating vapours.

Explosion hazard : May form flammable/explosive vapour-air mixture. Heavier than air, vapours may travel long

distances along ground, ignite and flash back to source.

5.3. Advice for firefighters

Firefighting instructions : Move containers away from the fire area if this can be done without risk. Cool closed

containers exposed to fire with water spray. Prevent runoff from entering water courses,

sewers and basements.

Protection during firefighting : Keep upwind of fire. Wear full fire fighting turn-out gear (full Bunker gear) and respiratory

protection (SCBA).

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General measures : Use personal protection recommended in Section 8. Isolate the hazard area and deny entry

to unnecessary and unprotected personnel. Use special care to avoid static electric

charges. Remove all sources of ignition.

6.1.1. For non-emergency personnel

No additional information available

6.1.2. For emergency responders

No additional information available

10/26/2023 (Revision date) EN (English) 3/13

Safety Data Sheet

For containment

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

6.2. Environmental precautions

Prevent entry to sewers and public waters.

6.3. Methods and material for containment and cleaning up

 Remove ignition sources. Stop leak if safe to do so. Absorb and/or contain spill with inert material (sand, vermiculite or other appropriate material), then place in suitable container.
 Do not flush into surface water or sewer system. Wear recommended personal protective

equipment.

Methods for cleaning up : Sweep or shovel spills into appropriate container for disposal. Provide ventilation.

6.4. Reference to other sections

For further information refer to section 8: "Exposure controls/personal protection".

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Additional hazards when processed

: Handle empty containers with care because residual vapours are flammable.

Precautions for safe handling

: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Avoid breathing dust/fume/gas/mist/vapours/spray. Avoid contact with skin, eyes and clothing. Do not swallow. Handle and open container with care. When using do not eat, drink or smoke. Take precautionary measures against static discharge. Use only non-

sparking tools. Use only outdoors or in a well-ventilated area.

Hygiene measures

: Take off immediately all contaminated clothing and wash it before reuse. Wash hands,

forearms and face thoroughly after handling.

7.2. Conditions for safe storage, including any incompatibilities

Technical measures Storage conditions : Proper grounding procedures to avoid static electricity should be followed.

 Keep out of the reach of children. Keep container tightly closed. Do not store in unlabelled containers. Store in dry, well-ventilated area. Keep cool. Keep out of direct sunlight.
 Containers which are opened should be properly resealed and kept upright to prevent

leakage. Protect from physical damage. Store locked up.

7.3. Specific end use(s)

Not available.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

8.1.1 National occupational exposure and biological limit values

Acetonitrile (75-05-8)		
EU - Indicative Occupational Exposure Limit (IOEL)		
IOEL TWA	70 mg/m³	
IOEL TWA [ppm]	40 ppm	
Remark	Possibility of significant uptake through the skin	
United Kingdom - Occupational Exposure Limits		
WEL TWA (OEL TWA) [1]	68 mg/m³	
WEL TWA (OEL TWA) [2]	40 ppm	
WEL STEL (OEL STEL)	102 mg/m³	
WEL STEL (OEL STEL) [ppm]	60 ppm	

Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

Formic acid (64-18-6)		
EU - Indicative Occupational Exposure Limit (IOEL)		
IOEL TWA	9 mg/m³	
IOEL TWA [ppm]	5 ppm	
United Kingdom - Occupational Exposure Limits		
WEL TWA (OEL TWA) [1]	9.6 mg/m³	
WEL TWA (OEL TWA) [2]	5 ppm	
WEL STEL (OEL STEL)	28.8 mg/m³ (calculated)	
WEL STEL (OEL STEL) [ppm]	15 ppm (calculated)	

8.1.2. Recommended monitoring procedures

Monitoring methods	
Monitoring methods	Consult the relevant monitoring standards for the region.

8.1.3. Air contaminants formed

No additional information available

8.1.4. DNEL and PNEC

No additional information available

8.1.5. Control banding

No additional information available

8.2. Exposure controls

8.2.1. Appropriate engineering controls

Appropriate engineering controls:

Ensure good ventilation of the work station. Provide readily accessible eye wash stations and safety showers.

8.2.2. Personal protection equipment

8.2.2.1. Eye and face protection

Eye protection:

Safety eyewear complying with an approved standard such as the European Standard EN166 should be used when a risk assessment indicates this is necessary to avoid exposure to liquid splashes, mists or dusts.

8.2.2.2. Skin protection

Skin and body protection:

Chemical resistant apron. Flame retardant and anti-static material recommended.

Hand protection:

Chemical resistant gloves (according to European standard NF ISO 374-1 or equivalent)

8.2.2.3. Respiratory protection

Respiratory protection:

In case of insufficient ventilation, wear suitable respiratory equipment. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator. SDSs cannot provide detailed and complete respiratory protection guidelines. Selection of respiratory protection must be done by a qualified person who has assessed the work environment.

8.2.2.4. Thermal hazards

Thermal hazard protection:

Not required for normal conditions of use.

8.2.3. Environmental exposure controls

Environmental exposure controls:

Avoid release to the environment.

Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

Other information:

Handle in accordance with good industrial hygiene and safety procedures. Do not eat, drink or smoke when using this product.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Liquid
Appearance : Transparent.
Colour : No data available
Odour : Pungent.

Odour threshold : No data available

pH : 2.7

Relative evaporation rate (butylacetate=1) : No data available
Melting point : No data available
Freezing point : No data available
Boiling point : 81 °C (177.8 °F)
Flash point : 2 °C (35.6 °F)
Auto-ignition temperature : No data available
Decomposition temperature : No data available

Flammability (solid, gas) : Highly flammable liquid and vapour.

: No data available Vapour pressure Relative vapour density at 20°C : No data available Relative density : No data available Solubility : Water: 100 % Partition coefficient n-octanol/water : No data available : No data available Viscosity, kinematic Viscosity, dynamic : No data available Explosive properties : No data available Oxidising properties : No data available **Explosive limits** : No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

No dangerous reactions known under normal conditions of use.

10.2. Chemical stability

Stable under normal conditions. May form flammable/explosive vapour-air mixture.

10.3. Possibility of hazardous reactions

No dangerous reactions known under normal conditions of use.

10.4. Conditions to avoid

Heat. Sparks. Open flame. Incompatible materials. Sources of ignition. Direct sunlight.

10.5. Incompatible materials

Strong acids. Strong bases. Strong oxidizers. Halogenated compounds. Perchlorates. Sulfites.

10.6. Hazardous decomposition products

May include, and are not limited to: oxides of carbon. Nitrogen oxides. Hydrogen cyanide. Ammonia. Acetic acid. May release flammable gases.

Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

SECTION 11: Toxicological information

SECTION 11: Toxicological information		
11.1 Information on toxicological effects		
Acute toxicity (oral) : Acute toxicity (dermal) : Acute toxicity (inhalation) :	Harmful if swallowed. Not classified. Not classified.	
ZipChip Peptides Diluent		
ATE CLP (oral)	318.149 mg/kg bodyweight	
Acetonitrile (75-05-8)		
LD50 dermal rabbit	> 2000 mg/kg	
LC50 inhalation rat	26.8 mg/l/4h	
Formic acid (64-18-6)		
LD50 oral rat	1100 mg/kg	
LD50 dermal rat	> 2000 mg/kg bodyweight Animal: rat, Guideline: OECD Guideline 402 (Acute Dermal Toxicity)	
LC50 inhalation rat	7.85 mg/l/4h	
Unknown acute toxicity (CLP) - SDS : Skin corrosion/irritation :	1% of the mixture consists of ingredient(s) of unknown acute toxicity (Oral) 1% of the mixture consists of ingredient(s) of unknown acute toxicity (Dermal) 1% of the mixture consists of ingredient(s) of unknown acute toxicity (Inhalation (Vapours)) Not classified. pH: 2.7	
Additional information :	Based on available data, the classification criteria are not met.	
Serious eye damage/irritation :	Causes serious eye irritation. pH: 2.7	
Respiratory or skin sensitisation : Additional information : Germ cell mutagenicity : Additional information : Carcinogenicity : Additional information : Additional information : Additional information :	Not classified. Based on available data, the classification criteria are not met. Not classified. Based on available data, the classification criteria are not met. Not classified. Based on available data, the classification criteria are not met.	
Formic acid (64-18-6)		
NOAEL (chronic, oral, animal/male, 2 years)	400 mg/kg bodyweight Animal: mouse, Animal sex: male, Guideline: OECD Guideline 453 (Combined Chronic Toxicity / Carcinogenicity Studies), Remarks on results: other:Effect type: toxicity (migrated information)	
Reproductive toxicity : Additional information : STOT-single exposure : Additional information : STOT-repeated exposure : Additional information :	Not classified. Based on available data, the classification criteria are not met. Not classified. Based on available data, the classification criteria are not met. Not classified. Based on available data, the classification criteria are not met.	
Acetonitrile (75-05-8)		
NOAEC (inhalation, rat, gas, 90 days)	400 ppm Animal: rat, Guideline: other:	
Formic acid (64-18-6)		
LOAEL (oral, rat, 90 days)	2000 mg/kg bodyweight Animal: rat, Guideline: OECD Guideline 453 (Combined Chronic Toxicity / Carcinogenicity Studies)	
NOAEL (oral, rat, 90 days)	400 mg/kg bodyweight Animal: rat, Guideline: OECD Guideline 453 (Combined Chronic Toxicity / Carcinogenicity Studies)	
NOAEC (inhalation, rat, dust/mist/fume, 90 days)	0.244 mg/l air Animal: rat, Guideline: OECD Guideline 413 (Subchronic Inhalation Toxicity: 90-Day Study)	

Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

Aspiration hazard : Not classified.

Additional information : Based on available data, the classification criteria are not met.

Other information : Likely routes of exposure: ingestion, inhalation, skin and eye

SECTION 12: Ecological information

12.1. Toxicity

Ecology - general : May cause long-term adverse effects in the aquatic environment.

Unknown hazards to the aquatic environment (CLP) : Contains 1 % of components with unknown hazards to the aquatic environment

Hazardous to the aquatic environment, short-term : Not classified.

(acute)

Hazardous to the aquatic environment, long-term

(chronic)

: Not classified.

Acetonitrile (75-05-8)	
LC50 - Fish [1]	1600 – 1690 mg/l (Exposure time: 96 h - Species: Pimephales promelas [flow-through])
LC50 - Fish [2]	1000 mg/l (Exposure time: 96 h - Species: Pimephales promelas [static])
EC50 - Crustacea [1]	> 1000 mg/l Test organisms (species): Daphnia magna
EC50 72h - Algae [1]	3560 mg/l Test organisms (species): Phaeodactylum tricornutum
EC50 72h - Algae [2]	9696 mg/l Test organisms (species): Phaeodactylum tricornutum
LOEC (chronic)	> 960 mg/l Test organisms (species): Daphnia magna Duration: '21 d'
NOEC (chronic)	960 mg/l Test organisms (species): Daphnia magna Duration: '21 d'
NOEC chronic fish	102 mg/l Test organisms (species): Oryzias latipes Duration: '21 d'
Formic acid (64-18-6)	
LC50 - Fish [1]	130 mg/l Test organisms (species): Danio rerio (previous name: Brachydanio rerio)
EC50 - Crustacea [1]	120 mg/l (Exposure time: 48 h - Species: Daphnia magna)
EC50 - Crustacea [2]	138 – 165.6 mg/l (Exposure time: 48 h - Species: Daphnia magna [Static])
EC50 72h - Algae [1]	26.9 mg/l (Species: Desmodesmus subspicatus)
EC50 96h - Algae [1]	25 mg/l (Species: Desmodesmus subspicatus)
LOEC (chronic)	> 100 mg/l Test organisms (species): Daphnia magna Duration: '21 d'
NOEC (chronic)	≥ 100 mg/l Test organisms (species): Daphnia magna Duration: '21 d'

12.2. Persistence and degradability

ZipChip Peptides Diluent	
Persistence and degradability	Not established.

12.3. Bioaccumulative potential

ZipChip Peptides Diluent		
Bioaccumulative potential	Not established.	
Acetonitrile (75-05-8)		
Partition coefficient n-octanol/water -0.34		
Formic acid (64-18-6)		
BCF - Fish [1]	(0.22 dimensionless)	
Partition coefficient n-octanol/water	-1.9 (at 23 °C (at pH 5)	

Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

PBT : No vPvB : No

ZipChip Peptides Diluent

This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII

This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII

12.6. Other adverse effects

Additional information : No other effects known

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Product/Packaging disposal recommendations : Dispose of contents/container to hazardous or special waste collection point, in accordance

with local, regional, national and/or international regulation. The generation of waste should

be avoided or minimized wherever possible.

Additional information : Handle empty containers with care because residual vapours are flammable.

SECTION 14: Transport information

In accordance with ADR / IMDG / IATA

14.1 UN number

UN-No. (ADR) : UN 1648 UN-No. (IMDG) : UN 1648 UN-No. (IATA) : UN 1648

14.2. UN proper shipping name

Proper Shipping Name (ADR) : ACETONITRILE
Proper Shipping Name (IMDG) : ACETONITRILE
Proper Shipping Name (IATA) : Acetonitrile

14.3. Transport hazard class(es)

ADR

Transport hazard class(es) (ADR) : 3

Danger labels (ADR)



IMDG

Transport hazard class(es) (IMDG) : 3
Danger labels (IMDG) : 3



Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

IATA

Transport hazard class(es) (IATA) : 3 Danger labels (IATA) 3



14.4. Packing group

Packing group (ADR) : 11 Packing group (IMDG) : 11 Packing group (IATA) : 11

14.5. Environmental hazards

Dangerous for the environment : No Marine pollutant No

Other information No supplementary information available.

14.6. Special precautions for user

Special transport precautions : Do not handle until all safety precautions have been read and understood.

Overland transport

Limited quantities (ADR) 11

Orange plates

33 1648

EAC code : •2YE

Transport by sea

Limited quantities (IMDG) : 1L

Air transport

PCA limited quantity max net quantity (IATA) : 1L PCA max net quantity (IATA) 5L

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

REACH Annex XVII (Restriction List)

Contains no substance(s) listed on REACH Annex XVII (Restriction Conditions)

REACH Annex XIV (Authorisation List)

Contains no substance(s) listed on REACH Annex XIV (Authorisation List)

REACH Candidate List (SVHC)

Contains no REACH candidate substance.

PIC Regulation (Prior Informed Consent)

Contains no substance(s) listed on the PIC list (Regulation EU 649/2012 concerning the export and import of hazardous chemicals)

POP Regulation (Persistent Organic Pollutants)

Contains no substance(s) listed on the POP list (Regulation EU 2019/1021 on persistent organic pollutants)

Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

Ozone Regulation (1005/2009)

Contains no substance(s) listed on the Ozone Depletion list (Regulation EU 1005/2009 on substances that deplete the ozone layer)

Explosives Precursors Regulation (2019/1148)

Contains no substance(s) listed on the Explosives Precursors list (Regulation EU 2019/1148 on the marketing and use of explosives precursors)

Drug Precursors Regulation (273/2004)

Contains no substance(s) listed on the Drug Precursors list (Regulation EC 273/2004 on the manufacture and the placing on market of certain substances used in the illicit manufacture of narcotic drugs and psychotropic substances)

15.1.2. National regulations

United Kingdom

British National Regulations : Not determined

15.2. Chemical safety assessment

No chemical safety assessment has been carried out

SECTION 16: Other information

Indication of changes			
Section	Changed item	Change	Comments
14	Transport information	Modified	V1.1
SDS	Name	Modified	V1.1

Abbreviations and acronyms:

°C – Degrees Celsius

°F – Degrees Fahrenheit

ADR - European Agreement concerning the International Carriage of Dangerous Goods by Road.

ACGIH - American Conference of Governmental Industrial Hygienists

ATE - Acute Toxicity Estimate

BCF - Bioconcentration Factor

BEI – Biological Exposure Index

CAS - Chemical Abstracts Service

CLP - Regulation (EC) No 1272/2008 on the Classification, Labeling and Packaging of substances and mixtures.

CMR - Carcinogen, Mutagen, Reproductive toxin

cP - centipoise (unit of dynamic viscosity)

cSt - centistokes (unit of kinematic viscosity)

DNEL - Derived No-effect Level

DMEL - Derived Minimal Effect Level

EC50 – Half maximal effective concentration

ECHA - European Chemicals Agency

EC-No. - European Community number

EU - European Union

GHS - Globally Harmonized System of Classification and Labelling of Chemicals

h – Hours

IATA - International Air Transport Association

IC50 - Inhibition concentration

IDLH – Immediately Dangerous to Life or Health

IMDG – International Maritime Dangerous Goods

IOELV - Indicative Occupational Exposure Limit Value

KIFS - Swedish Chemicals Agency's (Keml's) Code of Statutes

kPa - kilopascal

Koc - Adsorption Coefficient

Kow - Octanol-Water Partition Coefficient

LC50 - Median Lethal Concentration

LD50 - Median Lethal Dose

LOAEL - Lowest Observed Adverse Effect level

Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

Abbreviations and acronyms:

mg/l - Milligram per liter

mg/kg - Milligram per kilogram

mg/m3 - Milligram per cubic meter

Min - Minutes

NIOSH - National Institute for Occupational Safety and Health

NOEC - No Observed Effect Concentration

NO(A)EL - No Observed (Adverse) Effect Level

N.O.S. - Not Otherwise Specified

OEL - Occupational Exposure Limit

PBT - Persistent, Bioaccumulative and Toxic

PCN - Poison Centre Notification

PNEC - Predicted No Effect Concentration

ppm - Parts per million

PVC - Polyvinyl chloride

REACH - Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006

RID - European Agreement concerning the International Carriage of Dangerous Goods by Rail

SDS - Safety Data Sheet

STEL - Short Term Exposure Limit

STOT - Specific Target Organ Toxicity

SVHC - Substance of Very High Concern (CMR, vPvB, PBT)

TDI - Tolerable Daily Intake

TLV - Threshold Limit Value

TWA - Time Weighted Average

UFI - Unique Formulation Identifier

UN - United Nations

vPvB - Very Persistent and Very Bioaccumulative

WEL - Workplace Exposure Limit

WGK - Wassergefahrdungklasse - German water quality classification

Data sources : REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE

COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and

amending Regulation (EC) No 1907/2006.

Other information : None

Prepared by : Nexreg Compliance Inc.

www.Nexreg.com



Full text of H- and EUH-statements:	
Acute Tox. 3 (Inhalation:vapour)	Acute toxicity (inhalation:vapour) Category 3
Acute Tox. 4 (Dermal)	Acute toxicity (dermal), Category 4
Acute Tox. 4 (Inhalation)	Acute toxicity (inhal.), Category 4
Acute Tox. 4 (Inhalation:vapour)	Acute toxicity (inhalation:vapour) Category 4
Acute Tox. 4 (Oral)	Acute toxicity (oral), Category 4
Eye Dam. 1	Serious eye damage/eye irritation, Category 1
Eye Irrit. 2	Serious eye damage/eye irritation, Category 2
Flam. Liq. 2	Flammable liquids, Category 2
Flam. Liq. 3	Flammable liquids, Category 3
H225	Highly flammable liquid and vapour.
H226	Flammable liquid and vapour.
H302	Harmful if swallowed.

Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

Full text of H- and EUH-statements:	
H312	Harmful in contact with skin.
H314	Causes severe skin burns and eye damage.
H315	Causes skin irritation.
H318	Causes serious eye damage.
H319	Causes serious eye irritation.
H331	Toxic if inhaled.
H332	Harmful if inhaled.
Skin Corr. 1A	Skin corrosion/irritation, Category 1, Sub-Category 1A
Skin Corr. 1B	Skin corrosion/irritation, Category 1, Sub-Category 1B
Skin Irrit. 2	Skin corrosion/irritation, Category 2

Classification and procedure used to derive the classification for mixtures according to Regulation (EC) 1272/2008 [CLP]:		
Flam. Liq. 2	H225	On basis of test data
Acute Tox. 4 (Oral)	H302	Calculation method
Eye Irrit. 2	H319	Calculation method

Safety Data Sheet (SDS), EU - Nexreg Annex II 2022

Disclaimer: We believe the statements, technical information and recommendations contained herein are reliable, but they are given without warranty or guarantee of any kind. The information contained in this document applies to this specific material as supplied. It may not be valid for this material if it is used in combination with any other materials. It is the user's responsibility to satisfy oneself as to the suitability and completeness of this information for the user's own particular use.