

Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law Issue date: 3/2/2023 Revision date: 10/26/2023 Version: 1.1

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Mixture

Product name : ZipChip Amino Acid Standard

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category : Professional use

Use of the substance/mixture : For research and development use only.

Function or use category : Laboratory chemicals

1.2.2. Uses advised against

Restrictions on use : Not for use in diagnostic procedures.

1.3. Details of the supplier of the safety data sheet

Manufacturer

908 Devices 645 Summer St 02210 Boston, MA USA

T 1 (857) 254 - 1500 908devices.com

1.4. Emergency telephone number

Emergency number : 1 (844) 908 - 4357

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Safety data sheet available on request.

Full text of hazard classes, H- and EUH-statements: see section 16

Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

EUH-statements : EUH210 - Safety data sheet available on request.

Unknown acute toxicity (CLP) - SDS : 2% of the mixture consists of ingredient(s) of unknown acute toxicity (Dermal)

2% of the mixture consists of ingredient(s) of unknown acute toxicity (Inhalation (Dust/Mist))

2.3. Other hazards

Contains no PBT/vPvB substances ≥ 0.1% assessed in accordance with REACH Annex XIII

SECTION 3: Composition/information on ingredients

3.1. Substances

Not applicable

Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

3.2. Mixtures

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
2,3-Butanediol, 1,4-dimercapto-, (R*,R*)-	CAS-No.: 3483-12-3 EC-No.: 222-468-7	< 2	Acute Tox. 4 (Oral), H302 (ATE=300 mg/kg bodyweight) Skin Irrit. 2, H315 Eye Dam. 1, H318

Full text of H- and FUH-statements: see section 16

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures after inhalation : If breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable

for breathing. Get medical advice/attention if you feel unwell.

First-aid measures after skin contact : If skin irritation occurs: Wash skin with plenty of water. Obtain medical attention if irritation

persists.

First-aid measures after eye contact : IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if

present and easy to do. Continue rinsing. If eye irritation persists: Get medical

advice/attention.

First-aid measures after ingestion : Do not induce vomiting without medical advice. Never give anything by mouth to an

unconscious person. Get medical advice/attention if you feel unwell.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/effects after inhalation : May cause irritation to the respiratory tract.

Symptoms/effects after skin contact : May cause skin irritation. Repeated exposure may cause skin dryness or cracking.

Symptoms/effects after eye contact : May cause eye irritation. Symptoms may include discomfort or pain, excess blinking and

tear production, with possible redness and swelling.

Symptoms/effects after ingestion : May be harmful if swallowed. May cause gastrointestinal irritation, nausea, vomiting and

diarrhea.

4.3. Indication of any immediate medical attention and special treatment needed

Symptoms may be delayed. In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Carbon dioxide (CO2). Extinguishing powder. For large fire: Water spray. Alcohol resistant

foam.

Unsuitable extinguishing media : Do not use a solid water stream as it may scatter and spread fire.

5.2. Special hazards arising from the substance or mixture

Fire hazard : Products of combustion may include, and are not limited to: oxides of carbon.

5.3. Advice for firefighters

Firefighting instructions : Cool closed containers exposed to fire with water spray. Move containers away from the fire

area if this can be done without risk.

Protection during firefighting : Keep upwind of fire. Wear full fire fighting turn-out gear (full Bunker gear) and respiratory

protection (SCBA).

Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General measures

: Use personal protection recommended in Section 8. Isolate the hazard area and deny entry to unnecessary and unprotected personnel.

6.1.1. For non-emergency personnel

No additional information available

6.1.2. For emergency responders

No additional information available

6.2. Environmental precautions

Prevent entry to sewers and public waters.

6.3. Methods and material for containment and cleaning up

For containment

Methods for cleaning up

: Contain spill, then place in a suitable container. Do not flush to sewer or allow to enter waterways. Use appropriate Personal Protective Equipment (PPE).

-

: Sweep or shovel spills into appropriate container for disposal. Provide ventilation.

6.4. Reference to other sections

For further information refer to section 8: "Exposure controls/personal protection".

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling

: Avoid breathing dust. Avoid contact with skin and eyes. Do not swallow. When using do not eat, drink or smoke. Handle and open container with care. Wear appropriate PPE (see Section 8).

Hygiene measures

: Wash contaminated clothing before reuse. Always wash hands after handling the product.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions

: Keep out of the reach of children. Keep container tightly closed. Do not store in unlabelled containers. Store in dry, well-ventilated area. Keep cool. Keep out of direct sunlight. Containers which are opened should be properly resealed and kept upright to prevent leakage. Protect from physical damage.

7.3. Specific end use(s)

Not available.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

8.1.1 National occupational exposure and biological limit values

No additional information available

8.1.2. Recommended monitoring procedures

Monitoring methods

Monitoring methods Consult the relevant monitoring standards for the region.

8.1.3. Air contaminants formed

No additional information available

8.1.4. DNEL and PNEC

Additional information : Not applicable

10/26/2023 (Revision date) EN (English) 3/10

Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

8.1.5. Control banding

No additional information available

8.2. Exposure controls

8.2.1. Appropriate engineering controls

Appropriate engineering controls:

Ensure good ventilation of the work station.

8.2.2. Personal protection equipment

8.2.2.1. Eye and face protection

Eye protection:

Safety glasses or goggles are recommended when using product.

8.2.2.2. Skin protection

Skin and body protection:

Wear suitable protective clothing

Hand protection:

Wear suitable gloves. Consult glove manufacturer's product information on material suitability and material thickness.

8.2.2.3. Respiratory protection

Respiratory protection:

In case of insufficient ventilation, wear suitable respiratory equipment. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator. SDSs cannot provide detailed and complete respiratory protection guidelines. Selection of respiratory protection must be done by a qualified person who has assessed the work environment.

8.2.2.4. Thermal hazards

Thermal hazard protection:

Not required for normal conditions of use.

8.2.3. Environmental exposure controls

Environmental exposure controls:

Avoid release to the environment.

Other information:

Handle in accordance with good industrial hygiene and safety procedures. Do not eat, drink or smoke when using this product.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state: SolidAppearance: Transparent.Colour: No data availableOdour: No data availableOdour threshold: No data available

pH : 7

Relative evaporation rate (butylacetate=1) : No data available Melting point : No data available Freezing point : No data available Boiling point : No data available Flash point : No data available : No data available

Auto-ignition temperature : Product is not self-igniting

Decomposition temperature : No data available Flammability (solid, gas) Not flammable Vapour pressure : No data available : No data available Relative vapour density at 20°C Relative density No data available Solubility · Water: 100 % Partition coefficient n-octanol/water : No data available : No data available Viscosity, kinematic

10/26/2023 (Revision date) EN (English) 4/10

Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

Viscosity, dynamic : No data available
Explosive properties : No data available
Oxidising properties : No data available
Explosive limits : No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

No dangerous reactions known under normal conditions of use.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No dangerous reactions known under normal conditions of use.

10.4. Conditions to avoid

Heat

10.5. Incompatible materials

None known.

10.6. Hazardous decomposition products

May include, and are not limited to: oxides of carbon.

SECTION 11: Toxicological information

11.1 Information on toxicological effects

Acute toxicity (oral) : Not classified.
Acute toxicity (dermal) : Not classified.
Acute toxicity (inhalation) : Not classified.

Additional information : Based on available data, the classification criteria are not met.

2,3-Butanediol, 1,4-dimercapto-, (R*,R*)- (3483-12-3)

LD50 oral rat

300 – 2000 mg/kg bodyweight Animal: rat, Animal sex: female, Guideline: OECD
Guideline 423 (Acute Oral toxicity - Acute Toxic Class Method), Guideline: EU Method B.1
tris (Acute Oral Toxicity - Acute Toxic Class Method)

Unknown acute toxicity (CLP) - SDS : 2% of the mixture consists of ingredient(s) of unknown acute toxicity (Dermal)

2% of the mixture consists of ingredient(s) of unknown acute toxicity (Inhalation (Dust/Mist))

Skin corrosion/irritation : Not classified.

pH: 7

Additional information : Based on available data, the classification criteria are not met.

Serious eye damage/irritation : Not classified.

pH: 7

Additional information : Based on available data, the classification criteria are not met.

Respiratory or skin sensitisation : Not classified.

Additional information : Based on available data, the classification criteria are not met.

Germ cell mutagenicity : Not classified

Additional information : Based on available data, the classification criteria are not met.

Carcinogenicity : Not classifie

Additional information : Based on available data, the classification criteria are not met.

Reproductive toxicity : Not classified.

Additional information : Based on available data, the classification criteria are not met.

10/26/2023 (Revision date) EN (English) 5/10

Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

STOT-single exposure : Not classified.

Additional information : Based on available data, the classification criteria are not met.

STOT-repeated exposure : Not classified

Additional information : Based on available data, the classification criteria are not met.

Aspiration hazard : Not classified.

Additional information : Based on available data, the classification criteria are not met.

Other information : Likely routes of exposure: ingestion, inhalation, skin and eye

SECTION 12: Ecological information

12.1. Toxicity

Ecology - general : May cause long-term adverse effects in the aquatic environment.

Hazardous to the aquatic environment, short-term

(acute)

: Not classified.

Hazardous to the aquatic environment, long-term

: Not classified.

(chronic)

2,3-Butanediol, 1,4-dimercapto-, (R*,R*)- (3483-12-3)		
EC50 - Crustacea [1]	34.8 mg/l Test organisms (species): Daphnia magna	
EC50 72h - Algae [1]	24.3 mg/l Test organisms (species): Pseudokirchneriella subcapitata (previous names: Raphidocelis subcapitata, Selenastrum capricornutum)	
EC50 72h - Algae [2]	8.66 mg/l Test organisms (species): Pseudokirchneriella subcapitata (previous names: Raphidocelis subcapitata, Selenastrum capricornutum)	

12.2. Persistence and degradability

ZipChip Amino Acid Standard	
Persistence and degradability	Not established.

12.3. Bioaccumulative potential

ZipChip Amino Acid Standard	
Bioaccumulative potential	Not established.

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

ZipChip Amino Acid Standard

This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII

This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII

12.6. Other adverse effects

Additional information : No other effects known

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Product/Packaging disposal recommendations

: Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation. The generation of waste should be avoided or minimized wherever possible. Recycle empty containers where allowed.

Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

SECTION 14: Transport information

In accordance with ADR / IMDG / IATA

14.1 UN number

UN-No. (ADR) : Not regulated UN-No. (IMDG) : Not regulated UN-No. (IATA) : Not regulated

14.2. UN proper shipping name

Proper Shipping Name (ADR) : Not regulated Proper Shipping Name (IMDG) : Not regulated Proper Shipping Name (IATA) : Not regulated

14.3. Transport hazard class(es)

ADR

Transport hazard class(es) (ADR) : Not regulated

IMDG

Transport hazard class(es) (IMDG) : Not regulated

IATA

Transport hazard class(es) (IATA) : Not regulated

14.4. Packing group

Packing group (ADR) : Not regulated Packing group (IMDG) : Not regulated Packing group (IATA) : Not regulated

14.5. Environmental hazards

Dangerous for the environment : No Marine pollutant : No

Other information : No supplementary information available.

14.6. Special precautions for user

Special transport precautions : Do not handle until all safety precautions have been read and understood.

Overland transport

Not regulated

Transport by sea

Not regulated

Air transport

Not regulated

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

REACH Annex XVII (Restriction List)

Contains no substance(s) listed on REACH Annex XVII (Restriction Conditions)

Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

REACH Annex XIV (Authorisation List)

Contains no substance(s) listed on REACH Annex XIV (Authorisation List)

REACH Candidate List (SVHC)

Contains no REACH candidate substance.

PIC Regulation (Prior Informed Consent)

Contains no substance(s) listed on the PIC list (Regulation EU 649/2012 concerning the export and import of hazardous chemicals)

POP Regulation (Persistent Organic Pollutants)

Contains no substance(s) listed on the POP list (Regulation EU 2019/1021 on persistent organic pollutants)

Ozone Regulation (1005/2009)

Contains no substance(s) listed on the Ozone Depletion list (Regulation EU 1005/2009 on substances that deplete the ozone layer)

Explosives Precursors Regulation (2019/1148)

Contains no substance(s) listed on the Explosives Precursors list (Regulation EU 2019/1148 on the marketing and use of explosives precursors)

Drug Precursors Regulation (273/2004)

Contains no substance(s) listed on the Drug Precursors list (Regulation EC 273/2004 on the manufacture and the placing on market of certain substances used in the illicit manufacture of narcotic drugs and psychotropic substances)

15.1.2. National regulations

United Kingdom

British National Regulations

: Not determined.

15.2. Chemical safety assessment

No chemical safety assessment has been carried out

SECTION 16: Other information

Indication of changes			
Section	Changed item	Change	Comments
SDS	Name	Modified	V1.1

Abbreviations and acronyms:

°C – Degrees Celsius

°F – Degrees Fahrenheit

ADR - European Agreement concerning the International Carriage of Dangerous Goods by Road.

ACGIH – American Conference of Governmental Industrial Hygienists

ATE - Acute Toxicity Estimate

BCF - Bioconcentration Factor

BEI – Biological Exposure Index

CAS – Chemical Abstracts Service

CLP - Regulation (EC) No 1272/2008 on the Classification, Labeling and Packaging of substances and mixtures.

CMR – Carcinogen, Mutagen, Reproductive toxin

cP - centipoise (unit of dynamic viscosity)

cSt – centistokes (unit of kinematic viscosity)

DNEL - Derived No-effect Level

DMEL - Derived Minimal Effect Level

EC50 - Half maximal effective concentration

ECHA - European Chemicals Agency

EC-No. - European Community number

EU – European Union

GHS - Globally Harmonized System of Classification and Labelling of Chemicals

h – Hours

IATA – International Air Transport Association

Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

Abbreviations and acronyms:

IC50 - Inhibition concentration

IDLH - Immediately Dangerous to Life or Health

IMDG - International Maritime Dangerous Goods

IOELV - Indicative Occupational Exposure Limit Value

KIFS - Swedish Chemicals Agency's (Keml's) Code of Statutes

kPa - kilopascal

Koc - Adsorption Coefficient

Kow - Octanol-Water Partition Coefficient

LC50 - Median Lethal Concentration

LD50 - Median Lethal Dose

LOAEL - Lowest Observed Adverse Effect level

mg/l - Milligram per liter

mg/kg - Milligram per kilogram

mg/m3 - Milligram per cubic meter

Min - Minutes

NIOSH - National Institute for Occupational Safety and Health

NOEC - No Observed Effect Concentration

NO(A)EL - No Observed (Adverse) Effect Level

N.O.S. - Not Otherwise Specified

OEL - Occupational Exposure Limit

PBT - Persistent, Bioaccumulative and Toxic

PCN - Poison Centre Notification

PNEC - Predicted No Effect Concentration

ppm - Parts per million

PVC - Polyvinyl chloride

REACH - Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006

RID - European Agreement concerning the International Carriage of Dangerous Goods by Rail

SDS - Safety Data Sheet

STEL - Short Term Exposure Limit

STOT - Specific Target Organ Toxicity

SVHC - Substance of Very High Concern (CMR, vPvB, PBT)

TDI - Tolerable Daily Intake

TLV - Threshold Limit Value

TWA - Time Weighted Average

UFI - Unique Formulation Identifier

UN - United Nations

vPvB - Very Persistent and Very Bioaccumulative

WEL - Workplace Exposure Limit

WGK - Wassergefahrdungklasse - German water quality classification

Data sources : REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE

COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and

amending Regulation (EC) No 1907/2006.

Other information : None.

Prepared by : Nexreg Compliance Inc.

www.Nexreg.com

Full text of H- and EUH-statements:		
Acute Tox. 4 (Oral)	Acute toxicity (oral), Category 4	
EUH210	Safety data sheet available on request.	
Eye Dam. 1	Serious eye damage/eye irritation, Category 1	
H302	Harmful if swallowed.	
H315	Causes skin irritation.	
H318	Causes serious eye damage.	
Skin Irrit. 2	Skin corrosion/irritation, Category 2	

Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

Classification and procedure used to derive the classification for mixtures according to Regulation (EC) 1272/2008 [CLP]:		
EUH210	EUH210	Expert judgement

Safety Data Sheet (SDS), EU - Nexreg Annex II 2022

Disclaimer: We believe the statements, technical information and recommendations contained herein are reliable, but they are given without warranty or guarantee of any kind. The information contained in this document applies to this specific material as supplied. It may not be valid for this material if it is used in combination with any other materials. It is the user's responsibility to satisfy oneself as to the suitability and completeness of this information for the user's own particular use.