

# Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law Issue date: 1/6/2023 Revision date: 10/26/2023 Version: 1.1

# SECTION 1: Identification of the substance/mixture and of the company/undertaking

# 1.1. Product identifier

Product form : Mixture
Product name : Rebel BGE

## 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### 1.2.1. Relevant identified uses

Main use category : Professional use

Use of the substance/mixture : For research and development use only.

Function or use category : Laboratory chemicals

1.2.2. Uses advised against

Restrictions on use : Not for use in diagnostic procedures.

# 1.3. Details of the supplier of the safety data sheet

#### Manufacturer

908 Devices 645 Summer St 02210 Boston, MA USA

T 1 (857) 254 - 1500 908devices.com

### 1.4. Emergency telephone number

Emergency number : 1 (844) 908 - 4357

# **SECTION 2: Hazards identification**

# 2.1. Classification of the substance or mixture

### Classification according to Regulation (EC) No. 1272/2008 [CLP]

 Flam. Liq. 2
 H225

 Skin Irrit. 2
 H315

 Eye Irrit. 2
 H319

 STOT SE 3
 H336

Full text of hazard classes, H- and EUH-statements: see section 16

# Adverse physicochemical, human health and environmental effects

No additional information available

### 2.2. Label elements

### Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP)





GHS02

GHS07

Signal word (CLP) : Danger

Contains : Isopropyl alcohol; Formic acid

Hazard statements (CLP) : H225 - Highly flammable liquid and vapour.

H315 - Causes skin irritation.H319 - Causes serious eye irritation.H336 - May cause drowsiness or dizziness.

Precautionary statements (CLP) : P210 - Keep away from heat, hot surfaces, sparks, open flames and other ignition sources.

No smoking.

# Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

P233 - Keep container tightly closed.

P261 - Avoid breathing dust/fume/gas/mist/vapours/spray.

P264 - Wash hands, forearms and face thoroughly after handling.

P280 - Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.

P403+P235 - Store in a well-ventilated place. Keep cool.

### 2.3. Other hazards

Other hazards which do not result in classification : Exposure may aggravate those with pre-existing eye, skin, or respiratory conditions.

Contains no PBT/vPvB substances ≥ 0.1% assessed in accordance with REACH Annex XIII

# **SECTION 3: Composition/information on ingredients**

## 3.1. Substances

Not applicable

### 3.2. Mixtures

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
Isopropyl alcohol substance with national workplace exposure limit(s) (GB)	CAS-No.: 67-63-0 EC-No.: 200-661-7 EC Index-No.: 603-117-00-0	20 – 50	Flam. Liq. 2, H225 Eye Irrit. 2, H319 STOT SE 3, H336
Formic acid substance with national workplace exposure limit(s) (GB); substance with a Community workplace exposure limit (Note B)	CAS-No.: 64-18-6 EC-No.: 200-579-1 EC Index-No.: 607-001-00-0	< 4	Flam. Liq. 3, H226 Acute Tox. 4 (Oral), H302 (ATE=1100 mg/kg bodyweight) Acute Tox. 3 (Inhalation:vapour), H331 (ATE=7.85 mg/l/4h) Skin Corr. 1A, H314 Eye Dam. 1, H318

Specific concentration limits:		
Name	Product identifier	Specific concentration limits
Formic acid	CAS-No.: 64-18-6 EC-No.: 200-579-1 EC Index-No.: 607-001-00-0	$(2 \le C < 10)$ Skin Irrit. 2, H315 $(2 \le C < 10)$ Eye Irrit. 2, H319 $(10 \le C < 90)$ Skin Corr. 1B, H314 $(90 \le C < 100)$ Skin Corr. 1A, H314

Note B - Some substances (acids, bases, etc.) are placed on the market in aqueous solutions at various concentrations and, therefore, these solutions require different classification and labelling since the hazards vary at different concentrations. In Part 3 entries with Note B have a general designation of the following type: '... %'. In this case the supplier must state the percentage concentration of the solution on the label. Unless otherwise stated, it is assumed that the percentage concentration is calculated on a weight/weight basis.

Full text of H- and EUH-statements: see section 16

# **SECTION 4: First aid measures**

## 4.1. Description of first aid measures

First-aid measures after inhalation

: If inhaled and if breathing is difficult, remove victim to fresh air and keep at rest in a position comfortable for breathing. Call a POISON CENTER/doctor if you feel unwell.

First-aid measures after skin contact

: If on skin (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. Wash clothing before re-using. If skin irritation occurs: Get medical advice/attention.

10/26/2023 (Revision date) EN (English) 2/13

# Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if First-aid measures after eye contact

present and easy to do. Continue rinsing. If eye irritation persists: Get medical

advice/attention.

First-aid measures after ingestion Do not induce vomiting without medical advice. Never give anything by mouth to an

unconscious person. Get medical advice/attention if you feel unwell.

### 4.2. Most important symptoms and effects, both acute and delayed

Symptoms/effects after inhalation : May cause irritation to the respiratory tract. May cause drowsiness or dizziness.

Causes skin irritation. Symptoms may include redness, drying, defatting and cracking of the Symptoms/effects after skin contact

Symptoms/effects after eye contact : Causes serious eye irritation. Symptoms may include discomfort or pain, excess blinking

and tear production, with marked redness and swelling of the conjunctiva.

Symptoms/effects after ingestion May be harmful if swallowed. May cause gastrointestinal irritation, nausea, vomiting and

diarrhea

### 4.3. Indication of any immediate medical attention and special treatment needed

Symptoms may be delayed. In case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

# **SECTION 5: Firefighting measures**

### 5.1. Extinguishing media

Suitable extinguishing media : Dry chemical powder. Alcohol-resistant foam. Carbon dioxide (CO2). Unsuitable extinguishing media : Do not use a solid water stream as it may scatter and spread fire.

# 5.2. Special hazards arising from the substance or mixture

: Highly flammable liquid and vapour. Products of combustion may include, and are not Fire hazard

limited to: oxides of carbon. Irritating vapours. Formaldehyde.

Explosion hazard May form flammable/explosive vapour-air mixture. Heavier than air, vapours may travel long

distances along ground, ignite and flash back to source.

### 5.3. Advice for firefighters

Firefighting instructions : Move containers away from the fire area if this can be done without risk. Cool closed

containers exposed to fire with water spray. Prevent runoff from entering water courses,

sewers and basements.

Protection during firefighting : Keep upwind of fire. Wear full fire fighting turn-out gear (full Bunker gear) and respiratory

protection (SCBA).

### **SECTION 6: Accidental release measures**

### 6.1. Personal precautions, protective equipment and emergency procedures

General measures

: Use personal protection recommended in Section 8. Isolate the hazard area and deny entry to unnecessary and unprotected personnel. Use special care to avoid static electric charges. Remove all sources of ignition.

### 6.1.1. For non-emergency personnel

No additional information available

#### 6.1.2. For emergency responders

No additional information available

### 6.2. Environmental precautions

Prevent entry to sewers and public waters.

10/26/2023 (Revision date) EN (English) 3/13

# Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

## 6.3. Methods and material for containment and cleaning up

For containment : Stop leak if safe to do so. Remove ignition sources. Absorb and/or contain spill with inert

material (sand, vermiculite or other appropriate material), then place in suitable container. Do not flush into surface water or sewer system. Wear recommended personal protective

equipment.

Methods for cleaning up : Sweep or shovel spills into appropriate container for disposal. Provide ventilation.

## 6.4. Reference to other sections

For further information refer to section 8: "Exposure controls/personal protection".

# **SECTION 7: Handling and storage**

#### 7.1. Precautions for safe handling

Additional hazards when processed : Handle empty containers with care because residual vapours are flammable.

Precautions for safe handling : Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No

smoking. Avoid contact with skin and eyes. Avoid breathing

dust/fume/gas/mist/vapours/spray. Do not swallow. Handle and open container with care. When using do not eat, drink or smoke. Take precautionary measures against static discharge. Use only non-sparking tools. Use only outdoors or in a well-ventilated area.

: Take off immediately all contaminated clothing and wash it before reuse. Wash hands,

forearms and face thoroughly after handling.

#### 7.2. Conditions for safe storage, including any incompatibilities

Technical measures : Proper grounding procedures to avoid static electricity should be followed.

Storage conditions : Keep out of the reach of children. Keep container tightly closed. Do not store in unlabelled

containers. Store in a dry, cool and well-ventilated place. Keep cool. Keep out of direct sunlight. Containers which are opened should be properly resealed and kept upright to

prevent leakage. Protect from physical damage.

# 7.3. Specific end use(s)

Not available.

Hygiene measures

### **SECTION 8: Exposure controls/personal protection**

### 8.1. Control parameters

# 8.1.1 National occupational exposure and biological limit values

Isopropyl alcohol (67-63-0)		
United Kingdom - Occupational Exposure Limits		
Local name	Propan-2-ol	
WEL TWA (OEL TWA) [1]	999 mg/m³	
WEL TWA (OEL TWA) [2]	400 ppm	
WEL STEL (OEL STEL)	1250 mg/m³	
WEL STEL (OEL STEL) [ppm]	500 ppm	
Regulatory reference	EH40/2005 (Fourth edition, 2020). HSE	
Formic acid (64-18-6)		
EU - Indicative Occupational Exposure Limit (IOEL)		
IOEL TWA	9 mg/m³	
IOEL TWA [ppm]	5 ppm	

# Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

Formic acid (64-18-6)		
United Kingdom - Occupational Exposure Limits		
WEL TWA (OEL TWA) [1] 9.6 mg/m³		
WEL TWA (OEL TWA) [2]	5 ppm	
WEL STEL (OEL STEL)	28.8 mg/m³ (calculated)	
WEL STEL (OEL STEL) [ppm]	15 ppm (calculated)	

#### 8.1.2. Recommended monitoring procedures

Monitoring methods	
Monitoring methods	Consult the relevant monitoring standards for the region.

#### 8.1.3. Air contaminants formed

No additional information available

#### 8.1.4. DNEL and PNEC

No additional information available

#### 8.1.5. Control banding

No additional information available

## 8.2. Exposure controls

### 8.2.1. Appropriate engineering controls

### Appropriate engineering controls:

Ensure good ventilation of the work station. Provide readily accessible eye wash stations and safety showers.

### 8.2.2. Personal protection equipment

# 8.2.2.1. Eye and face protection

## Eye protection:

Safety eyewear complying with an approved standard such as the European Standard EN166 should be used when a risk assessment indicates this is necessary to avoid exposure to dusts. Contact lenses should not be worn

### 8.2.2.2. Skin protection

#### Skin and body protection:

Chemical resistant apron. Flame retardant and anti-static material recommended.

# Hand protection:

Chemical resistant gloves (according to European standard NF ISO 374-1 or equivalent)

# 8.2.2.3. Respiratory protection

# Respiratory protection:

In case of insufficient ventilation, wear suitable respiratory equipment. Respirator selection must be based on known or anticipated exposure levels, the hazards of the product and the safe working limits of the selected respirator. SDSs cannot provide detailed and complete respiratory protection guidelines. Selection of respiratory protection must be done by a qualified person who has assessed the work environment.

# 8.2.2.4. Thermal hazards

### Thermal hazard protection:

Not required for normal conditions of use.

### 8.2.3. Environmental exposure controls

# **Environmental exposure controls:**

Avoid release to the environment.

### Other information:

Handle in accordance with good industrial hygiene and safety procedures. Do not eat, drink or smoke when using this product.

# Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

### **SECTION 9: Physical and chemical properties**

### 9.1. Information on basic physical and chemical properties

Physical state : Liquid

Appearance : Clear, colourless liquid.

Colour : Colourless.

Odour Slight; Alcohol-like. Odour threshold No data available No data available pΗ Relative evaporation rate (butylacetate=1) No data available Melting point No data available Freezing point : No data available Boiling point No data available Flash point : No data available Auto-ignition temperature : No data available

Flammability (solid, gas) : Highly flammable liquid and vapour.

No data available

Solubility : Water: completely soluble

Partition coefficient n-octanol/water : No data available
Viscosity, kinematic : No data available
Viscosity, dynamic : No data available
Explosive properties : No data available
Oxidising properties : No data available
Explosive limits : No data available

### 9.2. Other information

Decomposition temperature

No additional information available

# **SECTION 10: Stability and reactivity**

# 10.1. Reactivity

No dangerous reactions known under normal conditions of use.

### 10.2. Chemical stability

Stable under normal conditions. May form flammable/explosive vapour-air mixture.

# 10.3. Possibility of hazardous reactions

No dangerous reactions known under normal conditions of use.

#### 10.4. Conditions to avoid

Heat. Incompatible materials. Sources of ignition. Direct sunlight.

# 10.5. Incompatible materials

Strong oxidizing agents. Aluminium. Magnesium. This material may attack some forms of plastics, rubbers and coatings.

### 10.6. Hazardous decomposition products

May include, and are not limited to: oxides of carbon. May release flammable gases. Formaldehyde.

## **SECTION 11: Toxicological information**

# 11.1 Information on toxicological effects

Acute toxicity (oral) : Not classified.

Acute toxicity (dermal) : Not classified.

# Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

Acute toxicity (inhalation)	Not classified.
Isopropyl alcohol (67-63-0)	
LD50 oral rat	5840 mg/kg bodyweight Animal: rat, Guideline: OECD Guideline 401 (Acute Oral Toxicity)
LD50 dermal rabbit	4059 mg/kg
LC50 inhalation rat	> 10000 ppm (Exposure time: 6 h)
Formic acid (64-18-6)	
LD50 oral rat	1100 mg/kg
LD50 dermal rat	> 2000 mg/kg bodyweight Animal: rat, Guideline: OECD Guideline 402 (Acute Dermal Toxicity)
LC50 inhalation rat	7.85 mg/l/4h
Skin corrosion/irritation Serious eye damage/irritation Respiratory or skin sensitisation Additional information Germ cell mutagenicity Additional information Carcinogenicity Additional information	Causes skin irritation. Causes serious eye irritation. Not classified. Based on available data, the classification criteria are not met. Not classified. Based on available data, the classification criteria are not met. Not classified. Based on available data, the classification criteria are not met. Based on available data, the classification criteria are not met.
Isopropyl alcohol (67-63-0)	
IARC group	3 - Not classifiable
Formic acid (64-18-6)	
NOAEL (chronic, oral, animal/male, 2 years)	400 mg/kg bodyweight Animal: mouse, Animal sex: male, Guideline: OECD Guideline 453 (Combined Chronic Toxicity / Carcinogenicity Studies), Remarks on results: other:Effect type: toxicity (migrated information)
Reproductive toxicity Additional information STOT-single exposure	Not classified. Based on available data, the classification criteria are not met. May cause drowsiness or dizziness.
Isopropyl alcohol (67-63-0)	
STOT-single exposure	May cause drowsiness or dizziness.
STOT-repeated exposure Additional information	Not classified.  Based on available data, the classification criteria are not met.
Formic acid (64-18-6)	
LOAEL (oral, rat, 90 days)	2000 mg/kg bodyweight Animal: rat, Guideline: OECD Guideline 453 (Combined Chronic Toxicity / Carcinogenicity Studies)
NOAEL (oral, rat, 90 days)	400 mg/kg bodyweight Animal: rat, Guideline: OECD Guideline 453 (Combined Chronic Toxicity / Carcinogenicity Studies)
NOAEC (inhalation, rat, dust/mist/fume, 90 days)	0.244 mg/l air Animal: rat, Guideline: OECD Guideline 413 (Subchronic Inhalation Toxicity: 90-Day Study)
Aspiration hazard Additional information Other information	Not classified. Based on available data, the classification criteria are not met. Likely routes of exposure: ingestion, inhalation, skin and eye

# SECTION 12: Ecological information

# 12.1. Toxicity

Ecology - general : May cause long-term adverse effects in the aquatic environment.

Hazardous to the aquatic environment, short–term : Not classified.

(acute)

10/26/2023 (Revision date) EN (English) 7/13

# Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

Hazardous to the aquatic environment, long-term : Not classified.

(chronic)

Isopropyl alcohol (67-63-0)	
LC50 - Fish [1]	10000 mg/l Test organisms (species): Pimephales promelas
LC50 - Fish [2]	9640 mg/l Test organisms (species): Pimephales promelas
EC50 - Crustacea [1]	13299 mg/l (Exposure time: 48 h - Species: Daphnia magna)
EC50 72h - Algae [1]	> 1000 mg/l (Species: Desmodesmus subspicatus)
EC50 96h - Algae [1]	> 1000 mg/l (Species: Desmodesmus subspicatus)
Formic acid (64-18-6)	
LC50 - Fish [1]	130 mg/l Test organisms (species): Danio rerio (previous name: Brachydanio rerio)
EC50 - Crustacea [1]	120 mg/l (Exposure time: 48 h - Species: Daphnia magna)
EC50 - Crustacea [2]	138 – 165.6 mg/l (Exposure time: 48 h - Species: Daphnia magna [Static])
EC50 72h - Algae [1]	26.9 mg/l (Species: Desmodesmus subspicatus)
EC50 96h - Algae [1]	25 mg/l (Species: Desmodesmus subspicatus)
LOEC (chronic)	> 100 mg/l Test organisms (species): Daphnia magna Duration: '21 d'
NOEC (chronic)	≥ 100 mg/l Test organisms (species): Daphnia magna Duration: '21 d'

# 12.2. Persistence and degradability

Rebel BGE	
Persistence and degradability	Readily biodegradable.

## 12.3. Bioaccumulative potential

Rebel BGE		
Bioaccumulative potential Not expected to acummulate.		
Isopropyl alcohol (67-63-0)		
Partition coefficient n-octanol/water 0.05 (at 25 °C)		
Formic acid (64-18-6)		
BCF - Fish [1] (0.22 dimensionless)		
Partition coefficient n-octanol/water	-1.9 (at 23 °C (at pH 5)	

# 12.4. Mobility in soil

No additional information available

# 12.5. Results of PBT and vPvB assessment

PBT : No vPvB : No

# Rebel BGE

This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII

This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII

## 12.6. Other adverse effects

Additional information : No other effects known

10/26/2023 (Revision date) EN (English) 8/13

# Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

### **SECTION 13: Disposal considerations**

### 13.1. Waste treatment methods

Product/Packaging disposal recommendations : Dispose of contents/container to hazardous or special waste collection point, in accordance

with local, regional, national and/or international regulation. The generation of waste should

be avoided or minimized wherever possible.

Additional information : Handle empty containers with care because residual vapours are flammable.

# **SECTION 14: Transport information**

In accordance with ADR / IMDG / IATA

### 14.1 UN number

UN-No. (ADR) : UN 1219 UN-No. (IMDG) : UN 1219 UN-No. (IATA) : UN 1219

## 14.2. UN proper shipping name

Proper Shipping Name (ADR) : ISOPROPANOL SOLUTION
Proper Shipping Name (IMDG) : ISOPROPANOL SOLUTION
Proper Shipping Name (IATA) : Isopropanol Solution

### 14.3. Transport hazard class(es)

#### **ADR**

Transport hazard class(es) (ADR) : 3
Danger labels (ADR) : 3



### **IMDG**

Transport hazard class(es) (IMDG) : 3
Danger labels (IMDG) : 3



### IATA

Transport hazard class(es) (IATA) : 3
Danger labels (IATA) : 3



# 14.4. Packing group

Packing group (ADR) : II
Packing group (IMDG) : II
Packing group (IATA) : II

### 14.5. Environmental hazards

Dangerous for the environment : No Marine pollutant : No

Other information : No supplementary information available.

# Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

### 14.6. Special precautions for user

Special transport precautions : Do not handle until all safety precautions have been read and understood.

Overland transport

Limited quantities (ADR) : 11 Excepted quantities (ADR) E2

Orange plates

33 : •2YE

EAC code

#### Transport by sea

No data available

#### Air transport

No data available

### 14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

## **SECTION 15: Regulatory information**

## 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### 15.1.1. EU-Regulations

### **REACH Annex XVII (Restriction List)**

Contains no substance(s) listed on REACH Annex XVII (Restriction Conditions)

### **REACH Annex XIV (Authorisation List)**

Contains no substance(s) listed on REACH Annex XIV (Authorisation List)

# **REACH Candidate List (SVHC)**

Contains no REACH candidate substance.

#### **PIC Regulation (Prior Informed Consent)**

Contains no substance(s) listed on the PIC list (Regulation EU 649/2012 concerning the export and import of hazardous chemicals)

### **POP Regulation (Persistent Organic Pollutants)**

Contains no substance(s) listed on the POP list (Regulation EU 2019/1021 on persistent organic pollutants)

### Ozone Regulation (1005/2009)

Contains no substance(s) listed on the Ozone Depletion list (Regulation EU 1005/2009 on substances that deplete the ozone layer)

### **Explosives Precursors Regulation (2019/1148)**

Contains no substance(s) listed on the Explosives Precursors list (Regulation EU 2019/1148 on the marketing and use of explosives precursors)

### **Drug Precursors Regulation (273/2004)**

Contains no substance(s) listed on the Drug Precursors list (Regulation EC 273/2004 on the manufacture and the placing on market of certain substances used in the illicit manufacture of narcotic drugs and psychotropic substances)

## 15.1.2. National regulations

Not determined

## 15.2. Chemical safety assessment

No chemical safety assessment has been carried out

# Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

### **SECTION 16: Other information**

Indication of changes			
Section Changed item Change Comments			
14	Transport information	Modified	V1.1
SDS	Name	Modified	V1.1

### **Abbreviations and acronyms:**

°C – Degrees Celsius

°F – Degrees Fahrenheit

ADR - European Agreement concerning the International Carriage of Dangerous Goods by Road.

ACGIH - American Conference of Governmental Industrial Hygienists

ATE – Acute Toxicity Estimate

BCF - Bioconcentration Factor

BEI – Biological Exposure Index

CAS - Chemical Abstracts Service

CLP - Regulation (EC) No 1272/2008 on the Classification, Labeling and Packaging of substances and mixtures.

CMR - Carcinogen, Mutagen, Reproductive toxin

cP - centipoise (unit of dynamic viscosity)

cSt - centistokes (unit of kinematic viscosity)

DNEL - Derived No-effect Level

DMEL - Derived Minimal Effect Level

EC50 - Half maximal effective concentration

ECHA - European Chemicals Agency

EC-No. - European Community number

EU - European Union

GHS - Globally Harmonized System of Classification and Labelling of Chemicals

h - Hours

IATA - International Air Transport Association

IC50 - Inhibition concentration

IDLH - Immediately Dangerous to Life or Health

IMDG – International Maritime Dangerous Goods

IOELV - Indicative Occupational Exposure Limit Value

KIFS - Swedish Chemicals Agency's (Keml's) Code of Statutes

kPa – kilopascal

Koc – Adsorption Coefficient

Kow - Octanol-Water Partition Coefficient

LC50 - Median Lethal Concentration

LD50 - Median Lethal Dose

LOAEL - Lowest Observed Adverse Effect level

mg/l - Milligram per liter

mg/kg - Milligram per kilogram

mg/m3 - Milligram per cubic meter

Min - Minutes

NIOSH - National Institute for Occupational Safety and Health

NOEC - No Observed Effect Concentration

NO(A)EL - No Observed (Adverse) Effect Level

N.O.S. - Not Otherwise Specified

OEL - Occupational Exposure Limit

PBT - Persistent, Bioaccumulative and Toxic

PCN - Poison Centre Notification

PNEC - Predicted No Effect Concentration

ppm - Parts per million

PVC - Polyvinyl chloride

REACH - Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006

RID - European Agreement concerning the International Carriage of Dangerous Goods by Rail

SDS - Safety Data Sheet

STEL - Short Term Exposure Limit

STOT - Specific Target Organ Toxicity

# Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

## **Abbreviations and acronyms:**

SVHC – Substance of Very High Concern (CMR, vPvB, PBT)

TDI - Tolerable Daily Intake

TLV - Threshold Limit Value

TWA - Time Weighted Average

UFI - Unique Formulation Identifier

UN - United Nations

vPvB - Very Persistent and Very Bioaccumulative

WEL - Workplace Exposure Limit

WGK - Wassergefahrdungklasse - German water quality classification

Data sources : REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE

COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and

amending Regulation (EC) No 1907/2006.

Other information : None.

Prepared by : Nexreg Compliance Inc.

www.Nexreg.com



Full text of H- and EUH-statements:			
Acute Tox. 3 (Inhalation:vapour)	Acute toxicity (inhalation:vapour) Category 3		
Acute Tox. 4 (Oral)	Acute toxicity (oral), Category 4		
Eye Dam. 1	Serious eye damage/eye irritation, Category 1		
Eye Irrit. 2	Serious eye damage/eye irritation, Category 2		
Flam. Liq. 2	Flammable liquids, Category 2		
Flam. Liq. 3	Flammable liquids, Category 3		
H225	Highly flammable liquid and vapour.		
H226	Flammable liquid and vapour.		
H302	Harmful if swallowed.		
H314	Causes severe skin burns and eye damage.		
H315	Causes skin irritation.		
H318	Causes serious eye damage.		
H319	Causes serious eye irritation.		
H331	Toxic if inhaled.		
H336	May cause drowsiness or dizziness.		
Skin Corr. 1A	Skin corrosion/irritation, Category 1, Sub-Category 1A		
Skin Corr. 1B	Skin corrosion/irritation, Category 1, Sub-Category 1B		
Skin Irrit. 2	Skin corrosion/irritation, Category 2		
STOT SE 3	Specific target organ toxicity – Single exposure, Category 3, Narcosis		

Classification and procedure used to derive the classification for mixtures according to Regulation (EC) 1272/2008 [CLP]:		
Flam. Liq. 2	H225	Calculation method
Skin Irrit. 2	H315	Specific concentration limits
Eye Irrit. 2	H319	Specific concentration limits
STOT SE 3	H336	Calculation method

Safety Data Sheet (SDS), EU - Nexreg Annex II 2022

# Safety Data Sheet

According to REACH Regulation (EC) No 1907/2006, as retained and amended in UK law

Disclaimer: We believe the statements, technical information and recommendations contained herein are reliable, but they are given without warranty or guarantee of any kind. The information contained in this document applies to this specific material as supplied. It may not be valid for this material if it is used in combination with any other materials. It is the user's responsibility to satisfy oneself as to the suitability and completeness of this information for the user's own particular use.